

Antibody Purification (small scale)

Vivian Liu

Abstract

This protocol described a easy procedure to purify anti-sera in small scale batch. Protein A or protein G beads can be used depending on the immunoglobulin types.

Citation: Vivian Liu Antibody Purification (small scale). **protocols.io**

dx.doi.org/10.17504/protocols.io.dny5fv

Published: 21 Aug 2015

Protocol

Step 1.

Centrifuge anti-sera @ 4°C, 14,000g, 10min.

📌 NOTES

Vivian Liu 20 Aug 2015

This step is to remove particulates in the serum.

■ ANNOTATIONS

Doug Holbrook 04 Sep 2015

Some more information

Step 2.

(Optional) Filter through 0.2µm filter.

📌 NOTES

Vivian Liu 20 Aug 2015

This step is to remove lipids in the serum.

Step 3.

Apply protein A or protein G beads gently to Bio-Rad empty gravity column.

📌 NOTES

Vivian Liu 20 Aug 2015

anti-sera:beads slurry=1:1(v/v)

Step 4.

Set up the column on a vacuum manifold. Equilibrate column with 10 volumes of ice cold Binding Buffer.

📌 NOTES

Vivian Liu 20 Aug 2015

Use vacuum to pack beads tightly.

Vivian Liu 21 Aug 2015

Binding Buffer: 1xPBS

Step 5.

Dilute anti-sera with 2 volumes of Binding Buffer.

Step 6.

Let diluted anti-sera go through the column by gravity.

📌 NOTES

Vivian Liu 20 Aug 2015

No vacuum here.

Step 7.

Collect flow through and repeat step 6 once.

Step 8.

Wash column with 10 bed volumes of ice-cold Binding Buffer.

Step 9.

Elute antibody with 5 bed volumes of room temperature Elution Buffer.

📌 NOTES

Vivian Liu 21 Aug 2015

Elution Buffer:

100mM citric acid, PH3.0 for protein A beads;

100mM Glycine-HCL, PH2.7 for protein G beads.

Step 10.

Fill up each collecting tube with 200 µl 1M Tris-HCl (PH9.0).

Step 11.

Collect 1ml of flow through in each collecting tube. Save for antibody activity test.

📌 NOTES

Vivian Liu 21 Aug 2015

The PH9.0 Tris-HCl in the collecting tube will take the PH of flow through to around 7.

Step 12.

Re-equilibrate column with 10 bed volumes of Binding Buffer. Store Column at refrigerator for later use.