

Cultivation of *Plasmodium falciparum* in serum-free media with growth-promoting factor from bovine plasma (GFS)

Version 2

Hiroko Asahi

Abstract

Citation: Hiroko Asahi Cultivation of *Plasmodium falciparum* in serum-free media with growth-promoting factor from bovine plasma (GFS). **protocols.io**

dx.doi.org/10.17504/protocols.io.i36cgre

Published: 23 Jul 2017

Protocol

Step 1.

step 1

The composition of media:

1. Complete culture medium (GFSRPMI):

Mix the following constituents aseptically.

- (1) 900ml RPMI1640 (with L-Glutamine, with 25mM HEPES, w/o NaHCO_3)*
- (2) 100ml GFS: GF21 (Wako Pure Chemical Industries, Osaka, Japan)
- (3) 2 g (24 mM) sodium bicarbonate (Invitrogen Ltd.)
- (4) 25mg gentamycin (Sigma-Aldrich Corp., St. Louis, MO, USA)
- (5) 20 mg (150 μM) hypoxanthine (Sigma-Aldrich)**

*gibco ref. 22400 (Invitrogen Ltd., Carlsbad, CA, USA)

**1.5 ml aliquot of 10 mM solution is added to 100 ml culture medium when using the medium.

2. Basal medium: Mix the following constituents aseptically.

- (1) 1000ml RPMI1640 (with L-Glutamine, with 25mM HEPES, w/o NaHCO_3)

(2) 25mg gentamycin (Sigma-Aldrich)

Step 2.

Step 2

Start cultures of *Plasmodium falciparum*

1. Preserve erythrocytes (RBCs) in Alsever's solution* at 4 °C for 3– 30 days.
2. Wash RBCs by centrifugation for 10 min at 800g, admixed with RBCs infected with *Plasmodium falciparum* (PfRBCs), dispense into 24-well culture plates at a hematocrit of 2% (1 ml of suspension/well), and culture in a humidified atmosphere of 5%CO₂, 5% O₂, and 90% N₂ at 37°C.
3. Adjust the parasitemia to 0.1% (for subculture) or 0.3%(for growth tests) by adding uninfected RBCs, and the hematocrit to 2% by adding the appropriate volume of culture medium.

*Alsever's solution:

1. Solubilize the following constituents in 200 ml distilled water.

Autoclave the mixture and add to whole blood or RBCs at 1:1 to 1:4.

- (1) Glucose 4.1g
- (2) NaCl 0.84g
- (3) Na₃ citrate.2H₂O 1.6g
- (4) Citric acid 0.11g

Step 3.

step 3

Preparation of highly, synchronized cultures at the ring stage:

The cultures of *Plasmodium falciparum* are synchronized at the ring stage by three successive exposures to 5% D-sorbitol.

1. Expose asynchronous PfRBCs to 5% (w/v) D-sorbitol (Sigma-Aldrich) for 5 min.
2. After washing with basal medium by centrifugation for 6 min at 500g, adjust the parasitemia to

0.8%–1.2% by adding uninfected RBCs, dispense into 24-well culture plates at a hematocrit of 2%, and cultured for 41 h.

3. At 41-h culture, expose the PfRBCs to 5% D-sorbitol (the 2nd exposure), wash, and culture for further 5 h.
4. At the total 46-h culture, expose the PfRBCs to 5% D-sorbitol (the 3rd exposure).
5. After the 3rd sorbitol treatment, remove residual schizonts and cell debris by isopycnic density centrifugation for 10 min at 800g on 63% Percoll PLUS* (GE Healthcare Bio-Science Corp., Tokyo, Japan).

*3.5 ml 90% Percoll PLUS + 1.5 ml PBS

90% Percoll PLUS is prepared by admixing 100% Percoll PLUS (9 parts) with 1.5 M NaCl (10x physiological saline) (1 part).

Note: Periods for growth to the schizont stage are somewhat variable among strains. Culture periods useful for the protocol should be confirmed beforehand.