

# Fasciola spp. faecal sedimentation protocol for concentration of eggs and DNA isolation Version 2

Nichola Calvani

## Abstract

**Citation:** Nichola Calvani Fasciola spp. faecal sedimentation protocol for concentration of eggs and DNA isolation. protocols.io

dx.doi.org/10.17504/protocols.io.jw3cpgn

**Published:** 18 Sep 2017

## Protocol

Traditional sedimentation for Fasciola spp.

### Step 1.

Mix faecal samples (3 g and 6 g for sheep and cattle, respectively) with distilled water to form a homogenous solution.

Traditional sedimentation for Fasciola spp.

### Step 2.

Hose the solution through a 270 µm nylon sieve into a 250 ml conical measuring cylinder, top with distilled water and allow to sediment for three minutes.

Traditional sedimentation for Fasciola spp.

### Step 3.

After three minutes aspirate the supernatant and pour the sediment into a 100 ml conical measuring cylinder, rinse the 250 ml conical cylinder into the new cylinder and top with distilled water. Allow to sediment for a further three minutes.

Traditional sedimentation for Fasciola spp.


### Step 4.

Aspirate the supernatant and pour the remaining sediment into a 15 ml centrifuge tube, rinse the 100 ml conical cylinder into the 15 ml centrifuge tube and top with distilled water. Allow to sediment for a final three minutes.

Traditional sedimentation for Fasciola spp.

### Step 5.


Aspirate the supernatant, leaving 2 ml of sediment. To perform a faecal egg count proceed to step 6. To go straight to DNA isolation proceed to step 9.



## Faecal egg count (FEC) for EPG calculation

### Step 6.

To examine the sediment for fluke eggs, add 2 drops of methylene blue (1%), shake to mix and rinse into a 6.5×17×1 cm grid perspex tray. Additional distilled water can be added to allow for ease of counting. Examine under a stereomicroscope at 15× magnification.



## Faecal egg count (FEC) for EPG calculation

### Step 7.

All yellow-brown *Fasciola* eggs should be counted. Counts should be divided by 2 to calculate eggs per gram (EPG) for cattle. EPGs for sheep are as observed.

## Concentration for disruption and DNA isolation

### Step 8.

If samples were used to calculate EPG: Return individual samples to the 15 ml centrifuge tubes and centrifuged at 2500 g for 10 minutes to form a pellet.

If proceeding straight from step 5: Centrifuge samples at 2500 g for 10 minutes to form a pellet.

#### Concentration for disruption and DNA isolation

##### **Step 9.**

Manually remove the entire pellet from the 15 ml centrifuge tube using a combination of Pasteur pipettes and fine wooden applicator sticks and place into a pre-prepared bead-beating tube containing ceramic beads and lysis buffer (BioLine Isolate Fecal DNA Kit).

#### Concentration for disruption and DNA isolation

##### **Step 10.**

Disrupt the samples at 6.0 m/s for 40 seconds on a bench top homogeniser (FastPrep®-24 MP Biomedicals, Australia). Place on ice after disruption until continuing with DNA isolation (BioLine Isolate Fecal DNA Kit) and amplification.