# **NC64A virus purification**

# **David Dunigan and Irina Agarkova**

# **Abstract**

**Citation:** David Dunigan and Irina Agarkova NC64A virus purification. **protocols.io** 

dx.doi.org/10.17504/protocols.io.etebeje

Published: 13 Jun 2016

## **Guidelines**

# **Supplemental notes:**

- 1) An  $A_{260}$  may be determined on a UV spectrophotometer (usually a 1:100 dilution works well).
  - 1 A<sub>260</sub> unit of PBCV-1 routinely yields 1.5-2.5 X 10<sup>10</sup> PFU/ml of virus.
- 2) For critical work, a second purification through sucrose gradients or a set of iodixanol gradients may be necessary.

#### **Protocol**

#### Step 1.

Inoculate flasks with Chlorella and incubate at 25°C with continuous light and shaking until the cells are in the actively growing phase. Inoculate flasks with *Chlorella* NC64A in MBBM (or *Micractinium* Pbi in FES) and incubate at  $25^{\circ}$ C with continuous light and shaking until the cells are in the actively growing phase (about 1-2 X  $10^{7}$  cells/ml).

# Step 2.

Infect the flasks of chlorella with virus at an moi of 0.01 to 0.001.

## Step 3.

Incubate the flasks for 48-72 hours at 25°C with continuous light and shaking. This material is now termed "lysate".

© DURATION

12:00:00

# Step 4.

Centrifuge the lysate in the Sorvall GSA rotor in 250 ml bottles at 5,000 rpm (4,000 rcf), 5 min, 4°C.

© DURATION 00:05:00

#### Step 5.

Discard the pellets.

## Step 6.

Add Triton X-100 to the lysate supernatants for a final concentration of 1% (from a 10 or 20% stock).

## Step 7.

Centrifuge the lysate in the Beckman Type 19 225 mL ultracentrifuge rotor at 17,000 rpm (43,000 rcf), 50 min, at 4°C.

## **O** DURATION

00:50:00

#### NOTES

#### Irina Agarkova 31 Mar 2016

Alternative 1, centrifuge the lysate in Beckman Ti 50.2 rotors at 20,000 rpm (24,000 rcf), 60 min, 4°C. Alternative 2, centrifuge in a GSA type rotor >10,000 rpm (15k X rcf), 2 h, 4°C.

#### Step 8.

Discard the supernatants.

## Step 9.

Resuspend the virus pellets with a small volume of 50 mM Tris-HCl, pH 7.8 (approximately 1.0 mL per 100 mL of original lysate).

# Step 10.

Layer the virus suspension onto 100-400 mg/mL (10-40%) linear sucrose density gradients equilibrated with 50 mM Tris-HCl, pH 7.8, made up in Beckman SW28 rotor tubes (layer approximately 3-4 mL per gradient).

## **Step 11.**

Centrifuge the gradients in a Beckman SW28 rotor at 20,000 rpm (72,000 rcf<sub>max</sub>), 20 min, 4°C.

## **O DURATION**

00:20:00

#### NOTES

## Irina Agarkova 26 Apr 2016

The virus will be the major band about 1/2 to 2/3 deep in the gradient.

#### **Step 12.**

Remove the virus bands from the gradients with sterile bent needles and transfer to oak ridge 30 mL polypropylene centrifuge tubes.

#### **Step 13.**

Split the virus from 3 gradients between 2 tubes.

#### **Step 14.**

Slowly dilute the virus to the tube volume with 50 mM Tris-HCl, pH 7.8.

#### **Step 15.**

Centrifuge the tubes in Beckman Ti 50.2 rotor at 27,000 rpm (44,000  $rcf_{max}$ ), 3 hours, 4°C. Alternatively, dilute the virus from the gradients 10-fold with Tris buffer and centrifuge in the Type 19 rotor for 1 hour, 17,000 rpm, 4 C. A GSA type high speed rotor can be used for at 2 hours, 12,000 rpm.

# **O** DURATION

03:00:00

#### **P** NOTES

# Irina Agarkova 31 Mar 2016

Alternatively, dilute the virus from the gradients  $\sim$ 10-fold with Tris buffer and centrifuge in the Type 19 rotor for 1 hour, 17,000 rpm (43,200 rcf<sub>max</sub>), 4°C.

## **Step 16.**

Discard the supernatants.

# **Step 17.**

Gently wash the pellet and bottle with some 50 mM Tris, pH 7.8 buffer to wash residual sucrose away.

## **Step 18.**

Resuspend the virus pellets with a small volume of 50 mM Tris-HCl, pH 7.8.

## Step 19.

Store the virus at 4°C. Do not freeze.

# Step 20.

Filter sterilizition using a 0.45 µm cellulose acetate or other low protein binding filter is recommended.