



Nov 11,  
2019

## Fluorescence microscopy with the marine heterotrophic flagellate *Cafeteria roenbergensis*

Monica Berjon-Otero<sup>1</sup>, Sarah Duponchel<sup>1</sup>, Matthias Fischer<sup>1</sup>

<sup>1</sup>Max Planck Institute for Medical Research

1 Works for me dx.doi.org/10.17504/protocols.io.qjmduk6

Protist Research to Optimize Tools in Genetics (PROT-G)

Monica Berjon-Otero

### ABSTRACT

Fluorescence microscopy is essential to understand, between others, the cell organization and the life cycle of viruses that infect them. This protocol describes the steps necessary to fix and permeabilize *Cafeteria roenbergensis* cells allowing its staining with antibodies or chemical reactions as Click-IT reactions.

### MATERIALS

NAME	CATALOG #	VENDOR
1.5 mL Eppendorf tubes		
PBS		
Sterile glass slides		
Paraformaldehyde fixative: 4% paraformaldehyde in phosphate buffered saline (PBS)		
Centrifuge		
Coverslip		
0.5% Triton in PBS1x		
ELGA water		
0.45 µm 25 mm diameter Millipore HA filter	View	
0.2 µm 25 mm Whatman Anodisc filter	View	
Glycerol		
Multi-hole vacuum system		

### Remove bacteria from *Cafeteria roenbergensis* culture

- Follow the protocol: "[Lysozyme-based removal of bacteria from cultures of the marine heterotrophic flagellate \*Cafeteria roenbergensis\*](#)"



If the culture is to be treated with other reagents, for example L-Azidohomoalanine, dilute the cells without bacteria to the required density and apply the treatment.

- Aliquot  $5 \times 10^6$  *Cafeteria* cells per 1.5 ml microfuge tube for each microscopy sample

#### Fixation and permeabilization

- 3 Centrifuge for 5 min at 5,000 g, 20°C
- 4 Resuspend each pellet in 500 µL of 4% paraformaldehyde (PFA) in 1X PBS
- 5 Incubate for 20 min at 20°C
- 6 Centrifuge for 5 min at 7,000 g, 20°C



You can keep the sample overnight at 4 °C by resuspending the pellet with 250 µL 1X PBS. On the following day, centrifuge for 5 min at 5,000 g, 20°C and continue with the next step

- 7 Resuspend the pellets in 500 µL of 0.5% Triton in 1X PBS
- 8 Incubate for 15 min at 20°C

#### Filters preparation to Cafeteria immobilization

- 9 Clean the filter-support grids of a vacuum manifold system with pure water (e.g. ELGA, milliQ)
- 10 Add a drop of ELGA water on each grid of the vacuum system that you want to use and place a 0.45 µm Millipore support filter on top of it
- 11 Apply vacuum to remove the drop of water
- 12 Add another drop of ELGA water on top of the Millipore filter
- 13 Place a 0.2 µm, 25 mm Anodisc Whatman filter on top of the ELGA water drop
- 14 Apply vacuum with a pressure between 5-10 Hg to remove the water
- 15 If necessary, block or wash the filter with the required solution. Remove the solution by applying vacuum
- 16 Place the *Cafeteria* sample from step 8 on top of the Whatman filter

- 17 Apply vacuum with a pressure between 5-10 Hg to remove the permeabilization buffer and to immobilize the cells on the filter



Perform all the necessary blocking and washing steps following the same instructions

#### Sample staining

- 18 Place the filter with the *Cafeteria* sample on a piece of Parafilm. Add the staining solution of interest on top of the filter and incubate following the manufacturer's instructions
- 19 Place the filter on the vacuum system again
- 20 Apply vacuum with a pressure between 5-10 Hg to remove the remaining staining buffer
- 21 If required, wash the filter again with the buffer of interest
- 22 If desired, stain DNA by incubating with an appropriate reagent (DAPI, SYBR, Hoechst)
- 23 Apply vacuum with a pressure between 5-10 Hg to remove the remaining liquid

#### Sample mounting

- 24 Add 10  $\mu$ L of 70% glycerol on a glass slide and place the filter on top of the glycerol



The filter must be placed with the *Cafeteria* cells in the top!

- 25 Add 10  $\mu$ L of 70% glycerol on the coverslip and place it on top of the filter
- 26 Apply pressure on the sample during 20 min at room temperature
- 27 Seal the coverslip with nail polish



Samples can be kept at 4 °C or analyzed with the microscope immediately



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited