

# Transformation of competent E.coli cells with plasmid DNA

Lenny Teytelman

## Abstract

Transformation of heat-shock competent E. coli cells

**Citation:** Lenny Teytelman Transformation of competent E.coli cells with plasmid DNA. **protocols.io**

dx.doi.org/10.17504/protocols.io.ccfstm

**Published:** 08 Oct 2014

## Before start

For incubation on ice, make sure the tubes are standing in an ice-water mix, because without water, the cooling effect of ice is not reproducible due to the air between the ice fragments, especially if you have to incubate for a certain period of time.

## Protocol

### Step 1.

Thaw the appropriate amount of competent cells on ice.

### Step 2.

Pre-chill the required number of empty 1.5 ml microcentrifuge tubes.

### Step 3.

Pipet 50 µl aliquots of cells into the pre-chilled tubes.

### Step 4.

Add 5-10 µl of a ligation reaction mix or 5 ng of pure plasmid DNA to each tube. Mix gently!

### Step 5.

Incubate the tubes of ice for 30 min

 **DURATION**

00:30:00

### Step 6.

Heat shock the cells for 45 sec at 42°C

 **DURATION**

00:00:45

### Step 7.

Place the tubes immediately on ice for at least 2 min

 **DURATION**

00:02:00

### Step 8.

Add 800 µl of SOC medium to each tube and incubate for 1 hour at 37°C

 **DURATION**

01:00:00

**Step 9.**

Transfer the cultures to 1.5 ml microcentrifuge tubes and spin for 1 min at 6000 rpm.

 **DURATION**

00:01:00

**Step 10.**

Remove 800 µl of the supernatant and resuspend the pellet.

**Step 11.**

Plate out the suspension on a LB agar plate containing the appropriate antibiotic.

**Step 12.**

Incubate the plates overnight at 37°C.

**Warnings**

If you notice a significant drop in colony numbers after several transformations with plasmid DNA it's time to prepare fresh competent cells.