

HIFI-Barcode SOP - Assembling COI barcodes using high-throughput sequencing Version 3

Shanlin Liu,Chentao Yang,Chengran Zhou,Xin Zhou

Abstract

We developed an Illumina-based pipeline, HIFI-Barcode, to produce full-length COI barcodes from pooled PCR amplicons generated by individual specimens. Using indexed primer sets and high-throughput sequencing platform strategy, and optimized analysis pipeline, the analytical cost and chemistry cost will significantly be reduced. The new protocol includes DNA preparation, amplification, and data analysis pipeline.

Citation: Shanlin Liu,Chentao Yang,Chengran Zhou,Xin Zhou HIFI-Barcode SOP - Assembling COI barcodes using high-throughput sequencing. **protocols.io**

dx.doi.org/10.17504/protocols.io.k9icz4e

Published: 17 May 2018

Before start

Before starting this work, we recommend you read the paper we have published, "Shanlin Liu, Chentao Yang, Chengran Zhou, Xin Zhou; Filling reference gaps via assembling DNA barcodes using high-throughput sequencing—moving toward barcoding the world, *GigaScience*, Volume 6, Issue 12, 1 December 2017, Pages 1-8, <https://doi.org/10.1093/gigascience/gix104>", and if you have any question, just contact us.

Materials

10X PCR buffer by [Takara](#)

exTaq by [Takara](#)


Cap strips by [Applied Biosystems](#)

10mM dNTP mix by [Takara](#)

Protocol

Overview

Step 1.



DNA preparation

Step 2.

Individual Genomic DNA could be extracted using the Glass Fiber Plate method following manufacturer's protocol or other existing method.

http://ccdb.ca/site/wp-content/uploads/2016/09/CCDB_DNA_Extraction.pdf

DNA amplification

Step 3.

3.1 PCR reaction mixture

Consumables & Equipment for PCR amplification

1.10X PCR Buffer for exTaq Store at -20°C.

2.10mM dNTP mix. Store at -20°C

- 3.molecular grade ddH₂O
- 4.exTaq store at -20°C
- 5.Microplate(Eppendorf® plates)
- 6.Cap strips or sealing film
- 7.Eppendorf 5830R high speed centrifuge
- 8.primers.

Ninety-six pairs of different tags were added to both ends of a common COI barcode primer set (LCO1490 and HCO2198) with each tag containing 5 bps allowing for ≥ 2 bp differences from each other.

Table 1. Indexed Primer sequences

Primer 5'to3'	Primer 5'to3'
Rev001 AAAGCTAACTTCAGGGTGACCAAAAAATCA	For001 AAAGCGGTCAACAAATCATAAAGATATTGG
Rev002 AACAGTAACTTCAGGGTGACCAAAAAATCA	For002 AACAGGGTCAACAAATCATAAAGATATTGG
Rev003 AACCTTAACTTCAGGGTGACCAAAAAATCA	For003 AACCTGGTCAACAAATCATAAAGATATTGG
Rev004 AACTCTAACTTCAGGGTGACCAAAAAATCA	For004 AACTCGGTCAACAAATCATAAAGATATTGG
Rev005 AAGCATAAACTTCAGGGTGACCAAAAAATCA	For005 AAGCAGGTCAACAAATCATAAAGATATTGG
Rev006 AAGGTAACTTCAGGGTGACCAAAAAATCA	For006 AAGGTGGTCAACAAATCATAAAGATATTGG
Rev007 AAGTGTAACCTTCAGGGTGACCAAAAAATCA	For007 AAGTGGGTCAACAAATCATAAAGATATTGG
Rev008 AATGGTAACTTCAGGGTGACCAAAAAATCA	For008 AATGGGGTCAACAAATCATAAAGATATTGG
Rev009 ACACTTAACTTCAGGGTGACCAAAAAATCA	For009 ACACTGGTCAACAAATCATAAAGATATTGG
Rev010 ACAGATAAACTTCAGGGTGACCAAAAAATCA	For010 ACAGAGGTCAACAAATCATAAAGATATTGG
Rev011 ACCATTAACTTCAGGGTGACCAAAAAATCA	For011 ACCATGGTCAACAAATCATAAAGATATTGG
Rev012 ACCTATAAACTTCAGGGTGACCAAAAAATCA	For012 ACCTAGGTCAACAAATCATAAAGATATTGG
Rev013 ACGAATAAACTTCAGGGTGACCAAAAAATCA	For013 ACGAAGGTCAACAAATCATAAAGATATTGG
Rev014 ACGTTTAACTTCAGGGTGACCAAAAAATCA	For014 ACGTTGGTCAACAAATCATAAAGATATTGG
Rev015 ACTGTAACTTCAGGGTGACCAAAAAATCA	For015 ACTGTGGTCAACAAATCATAAAGATATTGG
Rev016 AGAACTAACTTCAGGGTGACCAAAAAATCA	For016 AGAACGGTCAACAAATCATAAAGATATTGG
Rev017 AGACATAAACTTCAGGGTGACCAAAAAATCA	For017 AGACAGGTCAACAAATCATAAAGATATTGG
Rev018 AGAGTTAACTTCAGGGTGACCAAAAAATCA	For018 AGAGTGGTCAACAAATCATAAAGATATTGG
Rev019 AGATGTAACTTCAGGGTGACCAAAAAATCA	For019 AGATGGGTCAACAAATCATAAAGATATTGG
Rev020 AGCTTTAACTTCAGGGTGACCAAAAAATCA	For020 AGCTTGGTCAACAAATCATAAAGATATTGG
Rev021 AGGATTAACTTCAGGGTGACCAAAAAATCA	For021 AGGATGGTCAACAAATCATAAAGATATTGG
Rev022 AGTAGTAACTTCAGGGTGACCAAAAAATCA	For022 AGTAGGGTCAACAAATCATAAAGATATTGG
Rev023 AGTCTTAACTTCAGGGTGACCAAAAAATCA	For023 AGTCTGGTCAACAAATCATAAAGATATTGG
Rev024 AGTGATAAACTTCAGGGTGACCAAAAAATCA	For024 AGTGAGGTCAACAAATCATAAAGATATTGG
Rev025 AGTTCTAACTTCAGGGTGACCAAAAAATCA	For025 AGTTCGGTCAACAAATCATAAAGATATTGG
Rev026 ATACCTAACTTCAGGGTGACCAAAAAATCA	For026 ATACCGGTCAACAAATCATAAAGATATTGG
Rev027 ATCACTAACTTCAGGGTGACCAAAAAATCA	For027 ATCACGGTCAACAAATCATAAAGATATTGG

Rev028	CAAAGTAACTTCAGGGTGACCAAAAAATCA	For028	CAAAGGGTCAACAAATCATAAAGATATTGG
Rev029	CAACTTAACTTCAGGGTGACCAAAAAATCA	For029	CAACTGGTCAACAAATCATAAAGATATTGG
Rev030	CAATCTAACTTCAGGGTGACCAAAAAATCA	For030	CAATCGGTCAACAAATCATAAAGATATTGG
Rev031	CAGAATAAACTTCAGGGTGACCAAAAAATCA	For031	CAGAAGGTCAACAAATCATAAAGATATTGG
Rev032	CATACTAACTTCAGGGTGACCAAAAAATCA	For032	CATACGGTCAACAAATCATAAAGATATTGG
Rev033	CATCATAACTTCAGGGTGACCAAAAAATCA	For033	CATCAGGTCAACAAATCATAAAGATATTGG
Rev034	CCAATTAACTTCAGGGTGACCAAAAAATCA	For034	CCAATGGTCAACAAATCATAAAGATATTGG
Rev035	CGATTTAACTTCAGGGTGACCAAAAAATCA	For035	CGATTGGTCAACAAATCATAAAGATATTGG
Rev036	CGTATTAACTTCAGGGTGACCAAAAAATCA	For036	CGTATGGTCAACAAATCATAAAGATATTGG
Rev037	CGTTATAAACTTCAGGGTGACCAAAAAATCA	For037	CGTTAGGTCAACAAATCATAAAGATATTGG
Rev038	CTAACTAACTTCAGGGTGACCAAAAAATCA	For038	CTAACGGTCAACAAATCATAAAGATATTGG
Rev039	CTACATAAACTTCAGGGTGACCAAAAAATCA	For039	CTACAGGTCAACAAATCATAAAGATATTGG
Rev040	CTATGTAACTTCAGGGTGACCAAAAAATCA	For040	CTATGGGTCAACAAATCATAAAGATATTGG
Rev041	CTCAATAAACTTCAGGGTGACCAAAAAATCA	For041	CTCAAGGTCAACAAATCATAAAGATATTGG
Rev042	CTGATTAACTTCAGGGTGACCAAAAAATCA	For042	CTGATGGTCAACAAATCATAAAGATATTGG
Rev043	CTGTATAAACTTCAGGGTGACCAAAAAATCA	For043	CTGTAGGTCAACAAATCATAAAGATATTGG
Rev044	CTTAGTAACTTCAGGGTGACCAAAAAATCA	For044	CTTAGGGTCAACAAATCATAAAGATATTGG
Rev045	CTTCTTAACTTCAGGGTGACCAAAAAATCA	For045	CTTCTGGTCAACAAATCATAAAGATATTGG
Rev046	GAAACTAACTTCAGGGTGACCAAAAAATCA	For046	GAAACGGTCAACAAATCATAAAGATATTGG
Rev047	GAACATAAACTTCAGGGTGACCAAAAAATCA	For047	GAACAGGTCAACAAATCATAAAGATATTGG
Rev048	GAATGTAACTTCAGGGTGACCAAAAAATCA	For048	GAATGGGTCAACAAATCATAAAGATATTGG
Rev049	GACTTTAACTTCAGGGTGACCAAAAAATCA	For049	GACTTGGTCAACAAATCATAAAGATATTGG
Rev050	GAGATTAACTTCAGGGTGACCAAAAAATCA	For050	GAGATGGTCAACAAATCATAAAGATATTGG
Rev051	GAGTATAAACTTCAGGGTGACCAAAAAATCA	For051	GAGTAGGTCAACAAATCATAAAGATATTGG
Rev052	GATAGTAACTTCAGGGTGACCAAAAAATCA	For052	GATAGGGTCAACAAATCATAAAGATATTGG
Rev053	GATCTTAACTTCAGGGTGACCAAAAAATCA	For053	GATCTGGTCAACAAATCATAAAGATATTGG
Rev054	GATGATAAACTTCAGGGTGACCAAAAAATCA	For054	GATGAGGTCAACAAATCATAAAGATATTGG
Rev055	GATTCTAACTTCAGGGTGACCAAAAAATCA	For055	GATTCGGTCAACAAATCATAAAGATATTGG
Rev056	GCAATAAACTTCAGGGTGACCAAAAAATCA	For056	GCAAAGGTCAACAAATCATAAAGATATTGG
Rev057	GCTATTAACTTCAGGGTGACCAAAAAATCA	For057	GCTATGGTCAACAAATCATAAAGATATTGG
Rev058	GCTTATAAACTTCAGGGTGACCAAAAAATCA	For058	GCTTAGGTCAACAAATCATAAAGATATTGG
Rev059	GGAATTAACTTCAGGGTGACCAAAAAATCA	For059	GGAATGGTCAACAAATCATAAAGATATTGG
Rev060	GGATATAAACTTCAGGGTGACCAAAAAATCA	For060	GGATAGGTCAACAAATCATAAAGATATTGG
Rev061	GGTTTTAACTTCAGGGTGACCAAAAAATCA	For061	GGTTTGGTCAACAAATCATAAAGATATTGG
Rev062	GTAGATAAACTTCAGGGTGACCAAAAAATCA	For062	GTAGAGGTCAACAAATCATAAAGATATTGG
Rev063	GTCATTAACTTCAGGGTGACCAAAAAATCA	For063	GTCATGGTCAACAAATCATAAAGATATTGG
Rev064	GTGAATAAACTTCAGGGTGACCAAAAAATCA	For064	GTGAAGGTCAACAAATCATAAAGATATTGG
Rev065	GTGTTTAACTTCAGGGTGACCAAAAAATCA	For065	GTGTTGGTCAACAAATCATAAAGATATTGG
Rev066	GTTACTAACTTCAGGGTGACCAAAAAATCA	For066	GTTACGGTCAACAAATCATAAAGATATTGG
Rev067	G TTCATAAACTTCAGGGTGACCAAAAAATCA	For067	G TTCAGGTCAACAAATCATAAAGATATTGG
Rev068	TAAGGTAACTTCAGGGTGACCAAAAAATCA	For068	TAAGGGGTCAACAAATCATAAAGATATTGG
Rev069	TACTGTAACTTCAGGGTGACCAAAAAATCA	For069	TACTGGGTCAACAAATCATAAAGATATTGG
Rev070	TAGGATAAACTTCAGGGTGACCAAAAAATCA	For070	TAGGAGGTCAACAAATCATAAAGATATTGG
Rev071	TAGTCTAACTTCAGGGTGACCAAAAAATCA	For071	TAGTCGGTCAACAAATCATAAAGATATTGG

Rev072	TATCGTAAACTTCAGGGTGACCAAAAAATCA	For072	TATCGGGTCAACAAATCATAAAGATATTGG
Rev073	TATGCTAAACTTCAGGGTGACCAAAAAATCA	For073	TATGCGGTCAACAAATCATAAAGATATTGG
Rev074	TCACATAAACTTCAGGGTGACCAAAAAATCA	For074	TCACAGGTCAACAAATCATAAAGATATTGG
Rev075	TCAGTTAAACTTCAGGGTGACCAAAAAATCA	For075	TCAGTGGTCAACAAATCATAAAGATATTGG
Rev076	TCATGTAAACTTCAGGGTGACCAAAAAATCA	For076	TCATGGGTCAACAAATCATAAAGATATTGG
Rev077	TCCAATAAACTTCAGGGTGACCAAAAAATCA	For077	TCCAAGGTCAACAAATCATAAAGATATTGG
Rev078	TCCTTTAAACTTCAGGGTGACCAAAAAATCA	For078	TCCTTGGTCAACAAATCATAAAGATATTGG
Rev079	TCGATTAAACTTCAGGGTGACCAAAAAATCA	For079	TCGATGGTCAACAAATCATAAAGATATTGG
Rev080	TCGTATAAACTTCAGGGTGACCAAAAAATCA	For080	TCGTAGGTCAACAAATCATAAAGATATTGG
Rev081	TCTCTTAAACTTCAGGGTGACCAAAAAATCA	For081	TCTCTGGTCAACAAATCATAAAGATATTGG
Rev082	TGAAGTAAACTTCAGGGTGACCAAAAAATCA	For082	TGAAGGGTCAACAAATCATAAAGATATTGG
Rev083	TGACTTAAACTTCAGGGTGACCAAAAAATCA	For083	TGACTGGTCAACAAATCATAAAGATATTGG
Rev084	TGAGATAAACTTCAGGGTGACCAAAAAATCA	For084	TGAGAGGTCAACAAATCATAAAGATATTGG
Rev085	TGCTATAAACTTCAGGGTGACCAAAAAATCA	For085	TGCTAGGTCAACAAATCATAAAGATATTGG
Rev086	TGGAATAAACTTCAGGGTGACCAAAAAATCA	For086	TGGAAGGTCAACAAATCATAAAGATATTGG
Rev087	TGTACTAAACTTCAGGGTGACCAAAAAATCA	For087	TGTACGGTCAACAAATCATAAAGATATTGG
Rev088	TGTCATAAACTTCAGGGTGACCAAAAAATCA	For088	TGTCAGGTCAACAAATCATAAAGATATTGG
Rev089	TGTGTTAAACTTCAGGGTGACCAAAAAATCA	For089	TGTGTGGTCAACAAATCATAAAGATATTGG
Rev090	TTACGTAAACTTCAGGGTGACCAAAAAATCA	For090	TTACGGGTCAACAAATCATAAAGATATTGG
Rev091	TTAGCTAAACTTCAGGGTGACCAAAAAATCA	For091	TTAGCGGTCAACAAATCATAAAGATATTGG
Rev092	TTCTCTAAACTTCAGGGTGACCAAAAAATCA	For092	TTCTCGGTCAACAAATCATAAAGATATTGG
Rev093	TTGACTAAACTTCAGGGTGACCAAAAAATCA	For093	TTGACGGTCAACAAATCATAAAGATATTGG
Rev094	TTGCATAAACTTCAGGGTGACCAAAAAATCA	For094	TTGCAGGTCAACAAATCATAAAGATATTGG
Rev095	TTGGTTAAACTTCAGGGTGACCAAAAAATCA	For095	TTGGTGGTCAACAAATCATAAAGATATTGG
Rev096	TTTCCTAAACTTCAGGGTGACCAAAAAATCA	For096	TTTCCGGTCAACAAATCATAAAGATATTGG

PCR reagents per 25 µl reaction:

#of reactions	1x	10x
ddH ₂ O	16.2ul	1620ul
10uM Primer forward*	1ul	100ul
10uM Primer reverse*	1ul	100ul
10X primer buffer	3ul	300ul
DNTP mix	2.5ul	250ul
exTaq	0.3ul	30ul
Total	24ul	2400ul
DNA template	1ul per l well	

General recommendations

- The use of filter tips is recommended for all PCR reagents to avoid contamination. Clean the bench top with alcohol before setting up reactions.

- Always use a sterile tip when removing exTaq polymerase and the other reagents from their tubes.
- Keep DNA templates (i.e. other PCR products) away from the PCR reagents while you are setting up the reaction mixes. Add DNA after all of the reagents have been returned to the freezer.
- Always include a sample without template as a negative control to check for contamination of the reagents. Include a positive control (a DNA sample that has amplified in the past) as well to test the effectiveness of the PCR reagents.

3.2 PCR thermocycle program

The amplification program includes a thermocycling profile of :

- 94°C for 1 min;
- 5 cycles of 94°C for 30 sec, 45°C for 40 sec, and an extension at 72°C for 1 min;
- followed by 35 cycles of 94°C for 30 sec, 51°C for 40 sec, and 72°C for 1 min;
- with a final extension at 72°C for 10 min, and finally holding at 12°C.

3.3 PCR product check

All amplicons could be visualized on a 1.2% 96 Agarose E-gel (Biowest Agarose).



REAGENTS

10X PCR buffer by [Takara](#)

10mM dNTP mix by [Takara](#)

exTaq by [Takara](#)

Sequencing

Step 4.

4.1 Hiseq or BGISEQ500

All PCR products from each plate should be pooled using 1 µl per sample resulting in two 96 µl mixtures.

PCR amplicons should be fragmented to construct a library of an insert-size of 250 bp and sequenced with a strategy of 150 PE.

4.2 Pacbio (optional)

- All PCR products from each plate should be pooled using 1 µl per sample resulting in two 96 µl mixtures.
- Sequencing using PacBio RS II.

BIOINFORMATICS DATA ANALYSIS

Step 5.

5.1 HIFI-Barcode-Hiseq assembly

HIFI-Barcode-Hiseq could be downloaded from <https://github.com/comery/HIFI-barcode-hiseq>.

Description

HIFIBarcode is used to produce full-length COI barcodes from pooled PCR amplicons generated by individual specimens.

Installation

Clone from github

```
$ git clone https://github.com/comery/HIFI-barcode-hiseq.git
```

```
$ tar -zxvf HIFIBarcode.v1.0.tar.gz
```

Source code

(1) Wrapper script: HIFIBarcode.v1.0.pl: the wrapper script to run other Perl scripts to do the work you choose.

(2) Main scripts: Seven Perl scripts in folder HIFIBarcode.V1.0/bin/, including:

- 1_split_extract.pl
- 2_uniqu_sort_cluster.Pro.pl
- 3_sep_extract_overlap.pl
- 4_cluster_fromend.pl
- 5_forgap_filling.pl
- 6_rename_kmer.pl
- 7_final.pl

(3) Published softwares:

VSEARCH v2.4.4. Rognes, Torbjørn, et al. 'VSEARCH: a versatile open source tool for metagenomics.' PeerJ 4 (2016): e258. COPE CMR v1.0.3. Liu, Binghang, et al. 'COPE: an accurate k-mer-based pair-end reads connection tool to facilitate genome assembly.' Bioinformatics 28.22 (2012): 2870-2874. SOAPBarcode. Liu, Shanlin, et al. 'SOAPBarcode: revealing arthropod biodiversity through assembly of Illumina shotgun sequences of PCR amplicons.' Methods in Ecology and Evolution 4.12 (2013): 1142-1150.

Pre-requisites

Examples

1: run wrapper script to get shell text and then sh runHIFIBarcode.sh to run HIFIBarcode

```
perl HIFIBarcode.v1.0.pl --fq1 test_1.fq --fq2 test_2.fq --index index_primer.txt --length 5 --cpunum 10  
--outdir test --outpre testout
```

2: run shell

```
sh test/runHIFIBarcode.sh
```

NOTE

The proposed method assemble short-read Illumina sequences based on k-mer sequence matches and such misassembly was not observed in their real data, but it's still possible. The pipeline also provides an additional note file with a suffix of "note.txt" with notes alerting users about the possibilities which have similar or same scores comparing to their best alternative.

LATEST RELEASE

Version 1.1 201709

5.2 HIFI-Barcode-Pacbio assembly

HIFI-Barcode-Pacbio could be downloaded from <https://github.com/comery/HIFI-barcode-pacbio>.

DESCRIPTION

HIFIBarcode is used to produce full-length COI barcodes from pooled PCR amplicons generated by individual specimens.

INSTALLATION

- Clone from github

```
$ git clone https://github.com/comery/HIFI-barcode-pacbio.git
```


- Or go to website <https://github.com/comery/HIFI-barcode-pacbio> and click 'Download ZIP'

Requirements

(1) software

- PicBio smrtanalysis download from <http://www.pacb.com/products-and-services/analytical-software/smrt-analysis/>

(2) programming language

- standard perl
- standard python(python2 is ok)

(3) perl module

- Bio::Perl(exactly Bio::Seq)

(4) main perl and python scripts in bin/

- primer_like_extract.pl
- cluster_count_passes_length.pl
- py | see more information to go to [https://github.com/PacificBiosciences/Bioinformatics-Training/wiki/Extracting-Reads-of-Insert-\(C-\)-number-of-passes](https://github.com/PacificBiosciences/Bioinformatics-Training/wiki/Extracting-Reads-of-Insert-(C-)-number-of-passes) or directly download from https://raw.githubusercontent.com/PacificBiosciences/Bioinformatics-Training/master/scripts/ccs_passes.py
- pl

DATA requirements:

(1) pacbio original H5 file input

- data/*.h5

(2) primers list

- experiment_data/primer.lst

for GGTCACAAATCATAAAGATATTGG
rev TAACTTCAGGGTGACCAAAAAATCA

(3) index(barcodes for identifying samples) list

- experiment_data/index.xls

001 AAAGC
002 AACAG

003 AACCT
004 AACTC
005 AAGCA
...

(4) samples_location.tab

- samples name and corresponding location in 96-cell plate

1 A01
2 B01
3 C01
4 D01
5 E01
...

Overview of steps

If you installed PacBio smrtanalysis, I suppose you get the file path of setup.sh, more about Pacbio Data : <http://www.pacb.com/wp-content/uploads/SMRT-Link-User-Guide-v4.0.0.pdf>

e.g: setup_path='/path/PicBio/smrtanalysis/current/etc/setup.sh'

step 1 extract CCS from h5 files

Input:

my_inputs.fofn (*my_inputs.fofn contains files list of Pacbio H5 file in 01.data/, like this-> .01.data/m170506_092957_42199_c101149142550000001823255607191735_s1_p0.1.bax.h5*)

Output:

log

data *.ccs.fasta *.ccs.fastq *.ccs.h5 reads_of_insert.fasta reads_of_insert.fastq slots.pickle

workflow

results

Run:

```
$ source /path/PicBio/smrtanalysis/current/etc/setup.sh
```

```
$ fofnToSmrtpipeInput.py my_inputs.fofn > my_inputs.xml
```

```
$ smrtpipe.py --params=settings.xml xml:input.xml
```

step 2 extract passes number from CCS h5 files

Input:

```
/data/*.ccs.h5
```

Output:

```
ccs_passes.lst
```

Run:

Note: Before run the scripts, please in sure that you have run source /path/PicBio/smrtanalysis/current/etc/setup.sh

```
$ python bin/ccs_passes.py data/*.ccs.h5 >ccs_passes.lst
```

step 3 filtering CCS by passes number (>15)

Input:

```
ccs_passes.lst
```

```
data/reads_of_insert.fasta
```

Output:

```
ccs_passes_15.fa
```

Run:

```
$ awk '$2>=15{print $1}' ccs_passes.lst >ccs_passes_15.lst
```

```
$ perl ./bin/fish_ccs.pl ccs_passes_15.lst data/reads_of_insert.fasta >ccs_passes_15.fa
```

step 4 assigning CCS to samples by index

Input:

```
experiment_data/primer.fa
```

```
experiment_data/index.xls
```

```
ccs_passes_15.fa
```

Output:

Note: 'outdir' name is up to you, here default value is '02.assignment'

02.assignment/

assign.log.txt

ccs.successfully_assigned.fa

check.ccs_passes_15.fa.log

Run:

```
$ perl ./bin/1.primer_like_extract.pl -p experiment_data/primer.fa -index experiment_data/index.xls -fa ccs_passes_15.fa -cm 2 -cg 1
```

step 5 clustering CCS of each sample to find best one**Input:**

ccs.successfully_assigned.fa

check.ccs_passes_15.fa.log

ccs_passes.lst

Output:

cluster.top1.fas

cluster.id.txt

cluster.all.fa

Run:

```
$ cd 02.assignment/
```

```
$ perl ../bin/2.cluster_count_passes_length.pl -ccs ccs.successfully_assigned.fa -pattern check.ccs_passes_15.fa.log -passes ../ccs_passes.lst
```

```
$ perl ../bin/change_name-location.pl cluster.top1.fas >hifi-barcode-pacbio.cluster.top1.fa
```

ALL DONE!

So, 'hifi-barcode-pacbio.cluster.top1.fa' is final result!