

untitled protocol

Yuan Yao

Abstract

Citation: Yuan Yao untitled protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.qaedsbe

Published: 22 May 2018

Materials



X-Gal, 100mg v3941 by Promega

Protocol

Step 1.

The two recombinant plasmids encoding the T25-X (or X-T25) and T18-Y (or Y-T18) hybrid proteins are transformed into competent BTH101 reporter cells.

Step 2.

Co-transform two plasmids (pKNT-topA and pUT-rstA) into BTH101 competent Cells.

Step 3.

Transformants are plated on selective plates and incubated at 30°C for 40 hours.