



Oct 16, 2018

Working

Manual dissection of the *Schistosoma mansoni* head and back end for transcriptomic analysis

Sarah Seyffert¹¹iGEM 2018 Heinrich-Heine-University[dx.doi.org/10.17504/protocols.io.uquevww](https://doi.org/10.17504/protocols.io.uquevww)

iGEM Duesseldorf 2018

 Sarah Seyffert 

ABSTRACT

M2 medium for co-culturing Cyanobacteria, such as *S. elongatus*, together with heterotrophic bacteria

Slightly modified from:

Weiss, Taylor L., Eric J. Young, and Daniel C. Ducat. "A synthetic, light-driven consortium of cyanobacteria and heterotrophic bacteria enables stable polyhydroxybutyrate production." *Metabolic engineering* 44 (2017): 236-245.

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

Important information

- 1 Always work under sterile conditions.
Recipe for 1L 1x M2 medium

Trace metal mix

- 2 1000x concentration:
 - H_3BO_3 ($2.86 \text{ g} \cdot \text{L}^{-1}$)
 - $\text{MnCl}_2 \cdot 4 \text{ H}_2\text{O}$ ($1.81 \text{ g} \cdot \text{L}^{-1}$)
 - $\text{ZnSO}_4 \cdot 7 \text{ H}_2\text{O}$ ($0.222 \text{ g} \cdot \text{L}^{-1}$)
 - $\text{Na}_2\text{MoO}_4 \cdot 2 \text{ H}_2\text{O}$ ($0.390 \text{ g} \cdot \text{L}^{-1}$)
 - $\text{Co}(\text{NO}_3)_2 \cdot 6 \text{ H}_2\text{O}$ ($0.049 \text{ g} \cdot \text{L}^{-1}$)

Ingredients

- 3 Add 10 ml 100x M2 (<https://www.protocols.io/view/100x-m2-medium-uqnevev>)
- 4 Add 10 ml 100x K_2HPO_4 (82g/l)
- 5 Add 1 ml 1000x trace metals

6 Add 1 ml 1000x Na_2CO_3 ($20 \text{ mg} \cdot \text{L}^{-1}$)

7 Add 1 g HEPES

8 Fill up to 1 l with MilliQ-water

9 The medium has to be autoclaved afterwards

After autoclaving

10 Add 1 ml of 1000x Fe(III) ammonium citrate ($6 \text{ mg} \cdot \text{L}^{-1}$)

11 Add 200 μl of 5000x $\text{CuSO}_4 \cdot 5 \text{ H}_2\text{O}$ ($395 \text{ ng} \cdot \text{mL}^{-1}$)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited