



Antibiotic sensitivity test for dinoflagellates and Chromera

Yoshihisa Hirakawa¹, Hirokazu Sakamoto¹, Takashi Shiratori¹, Elisabeth Hehenberger², Nick Irwin², Patrick Keeling²

¹University of Tsukuba, ²University of British Columbia

dx.doi.org/10.17504/protocols.io.vcke2uw

Protist Research to Optimize Tools in Genetics (PROT-G)



Yoshihisa Hirakawa

University of Tsukuba



PROTOCOL STATUS

Working

We use this protocol in our group and it is working

Algal culture

- 1 Seven dinoflagellate species (*Amphidinium carterae*, *Glenodinium foliaceum*, *Heterocapsa triquetra*, *Karlodinium veneficum*, *Prorocentrum minimum*, *Symbiodinium minutum*, and *Togula britannica*) and *Chromera velia* were grown in 75 mL plastic flasks with Daigo IMK (Nihon Pharmaceutical Co., Ltd.), f/2, or ESM medium. The culture condition was at 20°C, under white illumination (60 $\mu\text{mol photons/m}^2/\text{s}^{-1}$) on a 14:10 light:dark cycle.

Antibiotic treatment

- 2 Precultured cells were transferred to 24-well plates, and an antibiotic drug (puromycin or zeocin) was added with different concentrations (0, 50, 100, 300 $\mu\text{g/mL}$). The cell density was 2,000 to 40,000 cells/mL at the starting point. The culture plates were incubated under the same condition described above.

Cell monitoring

- 3 The number of living cells was counted under a microscope using hemacytometer during a week.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited