

THP1 cell line differentiation

Rabia Khan

Abstract

Citation: Rabia Khan THP1 cell line differentiation. protocols.io

dx.doi.org/10.17504/protocols.io.fhubj6w

Published: 20 Oct 2016

Protocol

Step 1.

Grow cells till confluent - cells like to grow in clumps

Media used - RPMI + 10% FBS + b-mercaptoethanol

Ref: http://www.ncbi.nlm.nih.gov/pubmed/17334670

Step 2.

Plate 50,000 c/ well in a 24 well plate in 500 ul

Add 10-50ng/ml of PMA for 72 hours

Step 3.

After 3 days - rest cells for 16hours

Step 4.

Gamma stimulate with 10-50U IFNg o/n - use for PAMP stimulation