



Jun 12,  
2019

K/2 Ian / K-ET

Ian Probert<sup>1</sup>

<sup>1</sup>CNRS-Sorbonne Université, Station Biologique, Place G. Tessier 29680 Roscoff FRANCE

Working

dx.doi.org/10.17504/protocols.io.s53eg8n

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France



## ABSTRACT

Modified from K medium by Ian Probert

## Reference

Keller, M.D., Selvin, R.C., Claus, W. & Guillard, R.R.L. 1987. Media for the culture of oceanic ultraphytoplankton. J. Phycol. 23:633–8.

- 1 ■ To 994 ml of seawater (pH 8.2, adjusted with NaOH) add:

Quantity	Compound	Stock solution (sterile)	Final conc. in K medium
0.5ml	NaNO <sub>3</sub>	48.9542g/litre H <sub>2</sub> O	288µM
0.5ml	NH <sub>4</sub> Cl *	0.535g/litre H <sub>2</sub> O	5µM
0.5ml	KH <sub>2</sub> PO <sub>4</sub>	4.8992g/litre H <sub>2</sub> O	18µM
0.5ml	FeEDTA solution	(see recipe below)	(see below)
0.5ml	Trace metal solution	(see recipe below)	(see below)
1.0ml	f/2 vitamin solution	(see recipe below)	(see below)

\* optional

### FeEDTA solution

- 2

- To 950ml distilled H<sub>2</sub>O add:

Quantity	Compound	Stock solution	Final conc. in K medium
4.3g	(Na)FeEDTA	-	5.85µM

- Make up to 1 litre with milliQ H<sub>2</sub>O, sterilize (filter 0.22µm) and store in fridge.

### Trace metal solution

- 3

- To 950ml distilled H<sub>2</sub>O add:

Quantity	Compound	Stock solution	Final conc. in K medium
37.22g	Na <sub>2</sub> EDTA.2H <sub>2</sub> O	-	50µM
1.0ml	CuSO <sub>4</sub> .5H <sub>2</sub> O	2.497g/litre H <sub>2</sub> O	0.005µM
1.0ml	Na <sub>2</sub> MoO <sub>4</sub> .2H <sub>2</sub> O	7.2585g/litre H <sub>2</sub> O	0.015µM

1.0ml	ZnSO <sub>4</sub> ·7H <sub>2</sub> O	23.0g/litre H <sub>2</sub> O	0.004μM
1.0ml	CoSO <sub>4</sub> ·7H <sub>2</sub> O	14.055g/litre H <sub>2</sub> O	0.025μM
1.0ml	MnCl <sub>2</sub> ·4H <sub>2</sub> O	178.11g/litre H <sub>2</sub> O	0.45μM
1.0ml	H <sub>2</sub> SeO <sub>3</sub>	1.29g/litre H <sub>2</sub> O	0.005μM
1.0ml	NiCl <sub>2</sub> ·6H <sub>2</sub> O	1.49 g/litre H <sub>2</sub> O	0.00314μM

- Make up to 1 litre with milliQ H<sub>2</sub>O, sterilize (filter 0.22μm) and store in fridge.

#### f/2 Vitamin solution

- 4
- To 950ml distilled H<sub>2</sub>O add

Quantity	Compound	Stock solution	Final conc. in K medium
1.0ml	Vit. B <sub>12</sub> (cyanocobalamin)	0.5g/litre H <sub>2</sub> O	0.37nM
1.0ml	Biotin	5.0mg/litre H <sub>2</sub> O	2.0nM
100.0mg	Thiamine HCl	-	0.3μM

- Make up to 1 litre with milliQ H<sub>2</sub>O, filter sterilize into plastic vials and store in freezer.

#### Sterilization of medium

- 5
- Optional: *Heat to 80°C for 2 hours and leave to cool – this should kill most organisms but should not chemically modify the medium too much*
  - Filter sterilize through 0.22μm filters (e.g. Millipore Steritop units) into sterile (autoclaved) polycarbonate bottles.

#### For K-ET

- 6
- Add 10-30 ml marine soil extract (ET)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited