

# Obtaining pure cyanophage stocks (plaque purification)

Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen

## Abstract

**Citation:** Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen Obtaining pure cyanophage stocks (plaque purification). **protocols.io**

dx.doi.org/10.17504/protocols.io.dqq5vv

**Published:** 08 Feb 2016

## Protocol

### Step 1.

Make a dilution series of the lysate (assume  $10^4$  to  $10^5$  PFU per mL in the plaque lysate)

#### NOTES

**Amy Chan** 02 Sep 2015

Use this to perform a second round of plaque assays to purify the phage.

### Step 2.

Repeat the plaque purification procedure 2 more times to ensure that the cyanophage isolated is clonal.

### Step 3.

Repeat the plaque purification procedure again.

### Step 4.

Finally, prepare a primary cyanophage stock using lysate from the final purification via method A "[Liquid Amplification](#)" **OR** method B "[Plate Amplification](#)"

### Step 5.

Method A: liquid amplification

#### PROTOCOL

#### . [Liquid Amplification](#)

CONTACT: [Amy Chan](#)

### Step 5.1.

Add some of the lysate to target host in liquid culture.

### Step 5.2.

After the culture has lysed, remove cell debris via centrifugation.

### Step 5.3.

Filter sterilize the stock.

### Step 5.4.

Store at 4°C until further analysis.

### Step 6.

Method B: plate amplification

## PROTOCOL

### . [Plate Amplification](#)

CONTACT: [Amy Chan](#)

#### **Step 6.1.**

Prepare plaque assays with a dilution series of lysate from the final purification.

#### NOTES

**Amy Chan** 03 Sep 2015

Plates with confluent lysis of the host lawn (typically ca.  $10^4$  PFUs) can then be used to obtain cyanophage stocks by elution of phages from the plates.

#### **Step 6.2.**

Add 5 mL sterile seawater to the plate.

#### **Step 6.3.**

Scrape off the top agar layer into the seawater.

#### **Step 6.4.**

Leave at 4°C overnight.

#### DURATION

18:00:00

#### **Step 6.5.**

Remove agar and cell debris by centrifugation.

#### **Step 6.6.**

Filter sterilize the stock.

#### **Step 6.7.**

Store at 4°C until further analysis.

#### **Step 7.**

Titer the final stock via plaque assay.

#### **Step 8.**

Cyanophage stocks stored at 4°C in the dark are stable for at least a year.