

# DNA Extraction Procedure Using SDS Version 2

Jiang Hu

## Abstract

**Citation:** Jiang Hu DNA Extraction Procedure Using SDS. [protocols.io](https://www.protocols.io)

[dx.doi.org/10.17504/protocols.io.jg4cjyw](https://dx.doi.org/10.17504/protocols.io.jg4cjyw)

**Published:** 22 Aug 2017

## Materials

- Ethyl alcohol, Pure 200 proof, for molecular biology E7023 by [Sigma Aldrich](#)
- ✓ Liquid Nitrogen by Contributed by users
- 2-Mercaptoethanol by [Sigma Aldrich](#)
- Buffer EB 19086 by [Qiagen](#)
- Chloroform:Isoamyl alcohol 24:1 C0549 by [Sigma Aldrich](#)
- 2 x 0.5ml LongLife™ Proteinase K [5mg/ml] 786-038 by [G-Biosciences](#)
- 2 x 0.5ml LongLife™ RNase [10U/?l] 786-040 by [G-Biosciences](#)
- 1 Liter STE Buffer [10X] (100mM Tris.HCl (pH8.0), 10mM EDTA, 1M NaCl) 786-569 by [G-Biosciences](#)
- 100g SDS (Sodium dodecyl sulfate) DG092 by [G-Biosciences](#)
- 1kg Tris Base RC-106 by [G-Biosciences](#)
- ✓ isopropyl alcohol by Contributed by users

## Protocol

### Preparation

#### Step 1.

Cut tissue and grind by liquid nitrogen.

#### Step 2.

Transfer grinded tissue to 15 ml tube.

#### Step 3.

Add 6ml TEN, 700µl 20% SDS (56°C preheat), 150µl beta mercaptoethanol and 150 µl proteinase K, rapidly mix.

#### AMOUNT

6 ml Additional info: TEN

☐ AMOUNT

700 µl Additional info: 20% SDS

☐ AMOUNT

150 µl Additional info: beta mercaptoethanol

☐ AMOUNT

150 µl Additional info: proteinase K

## Incubation

### Step 4.

Incubate homogenate for 2 h at 56 °C , gently blending for every 5 10 min.

### Step 5.

Cool down to room temperature.

### Step 6.

Add an equal volume of Tris saturated phenol and mix.

### Step 7.

Centrifuge at room temperature (□16 °C) with 14000 RPM for 10 min, then save supernatant.

### Step 8.

Add chloroform and isoamyl alcohol (24:1) to the supernatant, then mix.

### Step 9.

Centrifuge with 14000 RPM for 10 min, save supernatant.

### Step 10.

Resuspend nuclei pellet with 0.8 X volume of frozen isopropyl alcohol.

### Step 11.

Wash with cold ethanol 75%. (1/2)

### Step 12.

Wash with cold ethanol 75%. (2/2)

### Step 13.

Let it dry.

### Step 14.

Add 200 µl TEN and 2 µl RNase (100 mg/ml).

☐ AMOUNT

200 µl Additional info: TEN

☐ AMOUNT

2 µl Additional info: RNase (100 mg/ml)

**Step 15.**

37 °C for 30 min, incubation period precipitation with Tip dolly, precipitate dissolve completely.

**Step 16.**

Add buffer (TEN and 200 µl 20% SDS) and 40 µl Protease K up to 2 ml.

☐ AMOUNT

200 µl Additional info: 20% SDS

☐ AMOUNT

40 µl Additional info: Protease K

**Step 17.**

Incubate for 30 min at 56 °C.

**Step 18.**

Add chloroform and isoamyl alcohol (24:1) to supernatant, then mix.

**Step 19.**

Centrifuge with 14000 RPM for 10 min, save supernatant.

**Step 20.**

Resuspend nuclei pellet with 0.8 X volume of frozen isopropyl alcohol.

**Step 21.**

Wash twice with cold ethanol 75%. (1/2)

**Step 22.**

Wash twice with cold ethanol 75%. (2/2)

**Step 23.**

Let it dry.

**Step 24.**

Add 300 µl EB (pH 8.0) to dissolve.

☐ AMOUNT

300 µl Additional info: EB