

Apoptosis

Lenka Zarybnicka, Ales Tichy

Abstract

Determination of apoptotic cell in lymphocyte population using flow-cytometry

Citation: Lenka Zarybnicka, Ales Tichy Apoptosis. **protocols.io**

dx.doi.org/10.17504/protocols.io.mjhc4j6

Published: 22 Jan 2018

Protocol

Step 1.

Take 2 ml of peripheral blood with 300 μ L of heparin at room temperature

Step 2.

Dilute 1 mL of Erythrocyte-Lysing Reagent (DAKO, Glostrup, Denmark) by 19 mL of deionized water

Step 3.

Add 2 mL of diluted Erythrocyte-Lysing Reagent to 100 μ L of peripheral blood and mix immediately

Step 4.

Incubate it for 10 minutes at room temperature

Step 5.

Centrifuge the sample (300 g, 5 mins, 24°C), then remove the supernatant

Step 6.

Wash the sample twice in 5 mL of cold phosphate-buffered saline

Step 7.

Dilute 1 mL of 10X concentrated Binding Buffer (BD Pharmingen™, USA) by 9 mL of deionized water

Step 8.

Resuspend cells of the sample in 1X Binding Buffer (BD Pharmingen™, USA) at a concentration of 5×10^6 cells/mL

Step 9.

Transfer 100 μ L of the solution (5×10^5 cells) to a 5 mL culture tube

Step 10.

Add 5 μ L of FITC Annexin V APC (BD Pharmingen™, USA) and gently vortex the cells

Step 11.

Incubate it for 10 min on ice in the dark

Step 12.

Wash the sample by 5 mL of cold 1X Binding buffer at 300g for 5 minutes at 5°C

Step 13.

Add 5 μ L of Propidium Iodide (250 μ g/mL) (BD Pharmingen™, USA)

Step 14.

Analyze by flow cytometry immediately.

