# DNA extraction of Nephila clavipes using the Wizard Genomic DNA Purification kit (Promega) Version 2

# Luiz Filipe Bartoleti, Luiz Filipe Bartoleti

# **Abstract**

Protocol used to extract whole genome DNA from Nephila clavipes (Araneae: Araneidae) specimens.

**Citation:** Luiz Filipe Bartoleti, Luiz Filipe Bartoleti DNA extraction of Nephila clavipes using the Wizard Genomic DNA Purification kit (Promega). **protocols.io** 

dx.doi.org/10.17504/protocols.io.knxcvfn

Published: 06 Nov 2017

## **Protocol**

## Step 1.

Macerate 1-8 legs (around 50 mg of tissue) using liquid nitrogen.

# Step 2.

Add 600 µL of Nuclei Lysis solution and homogenize.

## **■** AMOUNT

600 µl Additional info: Nuclei Lysis solution

## Step 3.

Take the microtubes to the dry bath at 65°C for one hour.

#### **↓** TEMPERATURE

65 °C Additional info:

## Step 4.

Add 200 µL of Protein Precipitation solution and homogeníze.

# **■** AMOUNT

200 µl Additional info: Protein Precipitation solution

# Step 5.

Centrifugate the microtubes at 13000 rpm for 10 minutes.

## Step 6.

Transfer the supernatant to a new microtube and add 600 µL of Isopropyl alcohol.

## **■** AMOUNT

600 µl Additional info: Isopropyl alcohol

# Step 7.

Keep the tubes at -20°C for at least one hour.

## **▮** TEMPERATURE

-20 °C Additional info:

## Step 8.

Centrifugate the microtubes at 13000 rpm for 6 minutes at 4°C.

## **▮** TEMPERATURE

4 °C Additional info:

## Step 9.

Discard the supernatant and add 600  $\mu L$  of absolut etanol.

## **■** AMOUNT

600 µl Additional info: Etanol

# Step 10.

Homogenize gently for one minute.

# **Step 11.**

Centrifugate the microtubes at 13000 rpm for 3 minutes at 4°C.

## **▮** TEMPERATURE

4 °C Additional info:

## **Step 12.**

Discard the supernatant.

# **Step 13.**

Let the samples dry.

# Step 14.

Rehydrate with 50 µL of Tris-EDTA buffer.

## **■** AMOUNT

50 μl Additional info: Tris-EDTA buffer

Published: 06 Nov 2017