



Jan 06,
2020

PI staining of Arabidopsis seedling

Magdalena Julkowska¹

¹King Abdullah University of Science and Technology

1 *Works for me* dx.doi.org/10.17504/protocols.io.pejdjcn

Salt Lab KAUST



Magdalena Julkowska
King Abdullah University of Science and Technology ⚡

transfer the seedlings from agar plate into a solution containing the nutrients as well as 5-10 µg/mL of propidium iodide. Incubate for 20 to 30 min.

1

Examine the red staining under (confocal) microscope

2

3



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited