

Radioactive in vitro transcription

Anna Behle

Abstract

Protocol for radioactive labelling of RNA using T7-Polymerase and α -[32 P]-UTP.

Citation: Anna Behle Radioactive in vitro transcription. **protocols.io**

dx.doi.org/10.17504/protocols.io.gudbws6

Published: 23 Mar 2017

Protocol

Buffers and reagents required

Step 1.

- 5x Ribomax buffer: 400 mM Hepes-KOH, pH 7.5; 60 mM MgCl_2 ; 10 mM spermidine; 200 mM DTT

Before start:

Step 2.

Generate a PCR product using primers containing the T7-promoter. Gel-extract or clean up using column purification.

Pipetting scheme

Step 3.

PCR Product containing the T7-Promoter (100 ng - 1 μg)

Ribomax buffer	6 μL
Pyrophosphatase	1 μL
RNase inhibitor	0.5 μL
DTT (1 M)	2 μL
ATP 10 mM	1.8 μL
CTP 10 mM	1.8 μL
GTP 10 mM	1.8 μL
α -[32 P]-UTP	3 μL

T7 RNA-Polymerase 1 μ L

H₂O *ad* 30 μ L

Incubate reaction for at least 2 hours at 37°C.

🕒 DURATION

02:00:00

Clean-up

Step 4.

Remove excessive radioactive nucleotides by cleaning the probe using a G-50 column (GE-Healthcare).