

Measuring spectral reflectance and transmittance (350-2500 nm) of large leaves using the Spectra Vista Corporation (SVC) DC-R/T Integrating Sphere Version 5

Etienne Laliberté, Raymond Soffer

Abstract

Here we describe the standardised protocol used by the [Canadian Airborne Biodiversity Observatory](#) (CABO) to measure leaf spectral reflectance and transmittance, using the [Spectra Vista Corporation](#) (SVC) [DC-R/T Integrating Sphere](#) fitted to a portable full-range [SVC HR-1024i](#) field spectroradiometer. This standard version of our protocol describes the common case where an individual leaf is large enough to entirely cover the reflectance or transmission port of the integrating sphere. Briefly, six mature, healthy-looking and sunlit leaves from a canopy plant are selected for measurements of adaxial reflectance and transmittance. Reflectance measurements are referenced to a calibrated Spectralon® disk and corrected for stray light. Our leaf spectroscopy protocol builds from [that](#) of the [Carnegie Airborne Observatory](#).

Citation: Etienne Laliberté, Raymond Soffer Measuring spectral reflectance and transmittance (350-2500 nm) of large leaves using the Spectra Vista Corporation (SVC) DC-R/T Integrating Sphere. **protocols.io**

dx.doi.org/10.17504/protocols.io.p8pdrvn

Published: 04 Jun 2018

Guidelines

Handling Spectralon®

- **Do not touch Spectralon®** (e.g. sphere interior, reference disks, plugs) with your fingers.
- **Do not use canned air** to remove dust on the Spectralon® disk; canned air contains chemicals that can alter Spectralon®'s optical properties.
- **Do not attempt to clean Spectralon®** in the field, other than **blowing surface dust only on the Spectralon® reference disk or sphere plugs** using the Canless Air Duster System; cleaning Spectralon® requires a special procedure that should only be done in the lab.

Equipment

- Spectra Vista Corporation [HR-1024i](#) full-range (350-2500 nm) field spectroradiometer
- Spectra Vista Corporation 3-inch Spectralon® DC-R/T [Sphere](#)
- Semi-rugged laptop or PDA running the SVC Scan software
- [Canless Air Duster System O₂ Hurricane](#) (**never use canned air**) to remove dust from the surface of the Spectralon® reference disk
- Plastic containers with lids to temporarily store leaf samples during measurements (optional)

Consumables

- Nitrile gloves for handling leaves

Before start

1. Consult the user manual of the spectroradiometer and the integrating sphere to set up the instrument.
2. The instrument should be set up in the shade, sheltered as much as possible from the elements (definitely not exposed to rain).
3. All canopy plants selected for measurements should have already been tagged, identified, and georeferenced before spectroscopy measurements start.
4. The spectroscopist should be positioned as close as possible to the sampled plants to minimise time from collection to measurement.
5. The spectroscopist should be in a comfortable position and have enough room around the instrument to spread leaf samples around without the risk of mixing up individual leaves during handling.
6. **Six mature, fully-developped, healthy-looking leaves from the sunlit (>3 h per day of direct sunlight) portion of the canopy are selected** for spectral measurements from the bulk leaf sample (often one of a few branches). Leaves should be collected from the uppermost surface of the branch (i.e. receiving the most direct sunlight),
7. Leaves should be stored in a sealed, clear plastic bag (breathe into the bag before sealing it) and **brought immediately to the spectroscopist for measurements.**

Protocol

Instrument Set-Up

Step 1.

Install the integrating sphere onto the spectroradiometer.

📌 NOTES

Etienne Laliberté 30 May 2018

Follow the SVC integrating sphere manual p. 6-16.

Instrument Set-Up

Step 2.

Power the spectroradiometer and **warm up for >15 min.**

📌 NOTES

Etienne Laliberté 26 Apr 2018

Record the **time of day** when the spectroradiometer was started.

Instrument Set-Up

Step 3.

Power the integrating sphere lamp and **warm up for >5 min.**

Software Set-Up

Step 4.

Ensure that the spectroradiometer is connected to the computer via Bluetooth.

📌 NOTES

Etienne Laliberté 30 May 2018

The default Bluetooth password for SVC HR-1024i spectroradiometers is **hr1024i**.

Software Set-Up

Step 5.

Launch SVC Scan.



Software Set-Up

Step 6.

Select **auto-save mode** using *File > Data Options > Auto Save*.

Software Set-Up

Step 7.

Select *File > New* to create a new working folder in your shared '**spectra**' Google Drive CABO folder.

📌 NOTES

Etienne Laliberté 30 May 2018

If you do not yet have a shared Google Drive folder for your project, contact the CABO data manager to create one (etienne.laliberte@umontreal.ca or jeremy.goimard@umontreal.ca).


Etienne Laliberté 30 May 2018

Your shared CABO Google Drive folder should be named *YYYYstarted-YourLastName-Degree-University*, e.g. 2018-Girard-MSc-UdeM.

Software Set-Up

Step 8.

Create a new folder named *YYYY-MM-DD-SiteID-SpectroradiometerSerialNumber* (without spaces) within that '**spectra**' folder, and keep the default file base name (= *grMMDDYY*).



NOTES


Etienne Laliberté 30 May 2018

The base file name is the prefix for all files that will be saved within that folder. In auto-save mode, each file will be added an incremental step, e.g., *gr053018_0000.sig*, *gr053018_0001.sig*, etc.

Software Set-Up

Step 9.

Go to *Control > Setup Instrument...* to connect the spectroradiometer.



NOTES

Etienne Laliberté 30 May 2018

Different instruments may be configured under different virtual COM ports.

Software Set-Up

Step 10.

1. Choose the *RAW DN* optic
2. Select 5 sec scan time

3. Select *Auto Integration*
4. Ensure the date and time are correct
5. Check *Integration Scales Raw DN Data*

Setup

Computer Port Control

☒ List Only Configured Ports ☐ List All Possible Ports

COM4

Disconnect Connect

HR-1024i Settings

Optic: RAW DN

Scan Timing

☒ Specify Total Scan Time (Sec): 5 ☐ Specify Coadds: Si: 4 Swir1: 32 Swir2: 32

Integration Time (mSec)

☒ Auto Integration Si: 50 Swir1: 20 Swir2: 6

Date: 06/04/18 Time: 11:18:30

Stored Scans: 0

Target Photo Acquisition: None

☒ Integration Scales RAW DN Data

Cancel OK

Software Set-Up

Step 11.

In *Control* > *Overlap / Matching Settings*:

1. Select *Preserve Overlapped Detector Data*
2. Select *Matching Type* > *None*

Overlap / Matching Settings

Detector Overlap Controls

- ☒ Preserve Overlapped Detector Data
- ☐ Remove Overlapped Detector Data

Si->Swir1 Transition Wvl: 990

Swir1->Swir2 Transition Wvl: 1900

Detector Matching Controls

Matching Type

- ☒ None
- ☐ Radiance
- ☐ Reflectance

Matching Region

Beginning Wvl: 980

Ending Wvl: 1000

Matching Factor Limits

Min: 0.80

Max: 1.20

NIR-SWIR Overlap Algorithm

- ☒ Use The NIR-SWIR Overlap Algorithm Within The Matching Region

Set Default Values Cancel OK

Review Protocol Summary Diagram

Step 12.

Review the [document](#) summarising the different sphere configurations (A–E), and the scans that need to be recorded in each configuration.

Configuration A: Reflectance Mode, Reference

Step 13.

Position the lamp over the sphere **primary light entrance port**.

SAFETY INFORMATION

The lamp can get very hot. Grab it by the slotted heat shield.

NOTES

Etienne Laliberté 24 Apr 2018

Make sure lamp is secured in locked position.

Configuration A: Reflectance Mode, Reference

Step 14.

Check lamp alignment.

NOTES

Etienne Laliberté 04 Jun 2018

Use a thin piece of paper at the exit of the reflectance sample port (empty port) to ensure the light beam under-fills and is centered in the reflectance port. **If it is not, then proceed to lamp alignment** as described in the SVC integrating sphere user manual, p. 29-32.

■ ANNOTATIONS

Etienne Laliberté 07 Jun 2018

enlever 'the' avant 'to'

Configuration A: Reflectance Mode, Reference

Step 15.

Screw the tethered light trap on the **reflectance port** sample holder.

🔗 NOTES

Etienne Laliberté 24 Apr 2018

The light trap can stay on the sample holder for the entire measurement session.

Configuration A: Reflectance Mode, Reference

Step 16.

Screw the untethered light trap on the **transmission port** sample holder.

🔗 NOTES

Etienne Laliberté 26 Apr 2018

The light trap can stay on the sample holder for all measurements made in reflectance mode (configurations A–C).

Configuration A: Reflectance Mode, Reference

Step 17.

Place the tethered calibrated Spectralon® reflectance standard over the **reflectance port**.

🔗 NOTES

Etienne Laliberté 26 Apr 2018

Place the standard over the reflectance port so that the light beam shines directly on its reflective surface (i.e facing into the sphere).

Configuration A: Reflectance Mode, Reference

Step 18.

Position leaf #1 over the **transmission port** with its adaxial (upper) surface facing into the sphere.

📌 NOTES

Etienne Laliberté 26 Apr 2018


Position the leaf so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein.

Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration A: Reflectance Mode, Reference

Step 19.

Collect a '**Reference Scan**' in this configuration.



📌 NOTES

Etienne Laliberté 25 Apr 2018

This corresponds to the **reference radiance** in reflectance mode (R_{ref}). The reference data will be automatically saved in all successive target scan files until a new 'Reference Scan' is made.

Configuration A: Reflectance Mode, Reference

Step 20.

Collect a '**Target Scan**' in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 04 Jun 2018

This reference 'target' scan, immediately collected after the reference, should give relative reflectance around 1 (+ noise) for all wavelengths. It is saved in a separate file to facilitate post-processing of spectra.

Configuration B: Reflectance Mode, Stray light

Step 21.

Carefully remove leaf #1 from the transmission port sample holder.

Configuration B: Reflectance Mode, Stray light

Step 22.

Remove the tethered calibrated Spectralon® reflectance standard from the reflectance port.

Configuration B: Reflectance Mode, Stray light


Step 23.

Place the tethered calibrated Spectralon® reflectance standard over the **transmission port** sample holder.

Configuration B: Reflectance Mode, Stray light

Step 24.

Collect a '**Target Scan**' in this configuration and **save the file**.



📌 NOTES

Etienne Laliberté 25 Apr 2018

This corresponds to the **stray light radiance** in reflectance mode (R_{str}).

Configuration C: Reflectance Mode, Target

Step 25.

Position leaf #1 over the **reflectance port** with its adaxial (upper) surface facing into the sphere.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reference radiance. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.


Etienne Laliberté 24 Apr 2018

The light trap should remain on the reflectance port sample holder.

Configuration C: Reflectance Mode, Target

Step 26.

Collect a '**Target Scan**' for leaf #1 in this configuration and **save the file**.



📌 NOTES

Etienne Laliberté 25 Apr 2018

This corresponds to the **target radiance** in reflectance mode for leaf #1 ($R_{tar,1}$).

Configuration C: Reflectance Mode, Target

Step 27.

Carefully replace leaf #1 by leaf #2.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reference radiance. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration C: Reflectance Mode, Target

Step 28.

Collect a '**Target Scan**' for leaf #2 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 24 Apr 2018

This corresponds to the **target radiance** in reflectance mode for leaf #2 ($R_{tar,2}$).

Configuration C: Reflectance Mode, Target

Step 29.

Carefully replace leaf #2 by leaf #3.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reference radiance. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration C: Reflectance Mode, Target

Step 30.

Collect a '**Target Scan**' for leaf #3 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 24 Apr 2018

This corresponds to the **target radiance** in reflectance mode for leaf #3 ($R_{tar,3}$).

Configuration C: Reflectance Mode, Target

Step 31.

Carefully replace leaf #3 by leaf #4.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reference radiance. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration C: Reflectance Mode, Target

Step 32.

Collect a '**Target Scan**' for leaf #4 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 24 Apr 2018

This corresponds to the **target radiance** in reflectance mode for leaf #4 ($R_{tar,4}$).

Configuration C: Reflectance Mode, Target

Step 33.

Carefully replace leaf #4 by leaf #5.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reference radiance. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration C: Reflectance Mode, Target

Step 34.

Collect a '**Target Scan**' for leaf #5 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 24 Apr 2018

This corresponds to the **target radiance** in reflectance mode for leaf #5 ($R_{tar,5}$).

Configuration C: Reflectance Mode, Target

Step 35.

Carefully replace leaf #5 by leaf #6.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reference radiance. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration C: Reflectance Mode, Target

Step 36.

Collect a '**Target Scan**' for leaf #6 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 24 Apr 2018

This corresponds to the **target radiance** in reflectance mode for leaf #6 ($R_{tar,6}$).

Configuration A: Reflectance Mode, Reference

Step 37.

Remove leaf #6 from the reflectance sample port holder.

Configuration A: Reflectance Mode, Reference

Step 38.

Place the tethered calibrated Spectralon® reflectance standard over the **reflectance port**.

📌 NOTES

Etienne Laliberté 26 Apr 2018

Place the standard over the reflectance port so that the light beam shines directly on its reflective surface (i.e facing into the sphere).

Configuration A: Reflectance Mode, Reference

Step 39.

Position leaf #1 over the **transmission port** with its adaxial (upper) surface facing into the sphere.

📌 NOTES


Etienne Laliberté 26 Apr 2018

Position the leaf so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein.
Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration A: Reflectance Mode, Reference

Step 40.

Collect a **'Target Scan'** in this configuration and **save the file**.



📌 NOTES

Etienne Laliberté 08 May 2018

This second reference radiance scan is only used to assess the stability of the system in reflectance mode.

Configuration D: Transmittance Mode, Reference

Step 41.

Remove the tethered calibrated Spectralon® reflectance standard from the sphere reflectance port.

Configuration D: Transmittance Mode, Reference

Step 42.

Remove leaf #1 from the transmission port sample holder.

Configuration D: Transmittance Mode, Reference

Step 43.

Position leaf #1 over the **reflectance port** so that its **abaxial** (lower) side is now facing the inside of the sphere.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces the inside of the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration D: Transmittance Mode, Reference

Step 44.

Remove the light trap from the transmission port sample holder.

Configuration D: Transmittance Mode, Reference

Step 45.

Position the lamp over the sphere **transmission port**.

⚠ SAFETY INFORMATION

The lamp can get very hot. Grab it by the slotted heat shield.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Make sure lamp is secured in locked position.

Configuration D: Transmittance Mode, Reference

Step 46.

Install the Spectralon® plug over the **primary light port**.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Ensure that the curved plug is placed the correct way to match the curvature of the sphere.


Configuration D: Transmittance Mode, Reference

Step 47.

Collect a '**Reference Scan**' in this configuration.

Transmittance mode

Reference radiance (T_{ref})



⊕ NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **reference radiance** in **transmittance** mode (T_{ref}).

Configuration D: Transmittance Mode, Reference

Step 48.

Collect a '**Target Scan**' in this configuration and **save the file**.

⊕ NOTES

Etienne Laliberté 04 Jun 2018

This reference 'target' scan, immediately collected after the reference, should give relative reflectance around 1 (+ noise) for all wavelengths. It is saved in a separate file to facilitate post-processing of spectra.

Configuration E: Transmittance Mode, Target

Step 49.

Carefully remove leaf #1 from the reflectance port sample holder.

⊕ NOTES

Etienne Laliberté 24 Apr 2018

The reflectance port should now be **empty** (but with **light trap on**).

Configuration E: Transmittance Mode, Target

Step 50.

Gently pull lamp and transmission sample holder away from the sphere.

SAFETY INFORMATION

The lamp can get very hot. Grab it by the slotted heat shield.

Configuration E: Transmittance Mode, Target

Step 51.

Place leaf #1 over the **transmission port** with its **abaxial** (lower) surface facing into the sphere.

NOTES

Etienne Laliberté 26 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces into the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration E: Transmittance Mode, Target

Step 52.

Release the transmission sample holder and move lamp back to its locked position.


Configuration E: Transmittance Mode, Target

Step 53.

Collect a '**Target Scan**' for leaf #1 in this configuration and **save the file**.

Transmittance mode

Target radiance, leaf i ($T_{tar,i}$)



⊕ NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf #1 ($T_{tar,1}$).

Configuration E: Transmittance Mode, Target

Step 54.

Carefully replace leaf #1 by leaf #2.

⊕ NOTES

Etienne Laliberté 26 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces into the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration E: Transmittance Mode, Target

Step 55.

Collect a '**Target Scan**' for leaf #2 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf #2 ($T_{\text{tar},2}$).

Configuration E: Transmittance Mode, Target

Step 56.

Carefully replace leaf #2 by leaf #3.

📌 NOTES

Etienne Laliberté 25 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces the inside of the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration E: Transmittance Mode, Target

Step 57.

Collect a '**Target Scan**' for leaf #3 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf #3 ($T_{\text{tar},3}$).

Configuration E: Transmittance Mode, Target

Step 58.

Carefully replace leaf #3 by leaf #4.

📌 NOTES

Etienne Laliberté 25 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces the inside of the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration E: Transmittance Mode, Target

Step 59.

Collect a '**Target Scan**' for leaf #4 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf #4 ($T_{\text{tar},4}$).

Configuration E: Transmittance Mode, Target

Step 60.

Carefully replace leaf #4 by leaf #5.

📌 NOTES

Etienne Laliberté 25 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces the inside of the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration E: Transmittance Mode, Target

Step 61.

Collect a '**Target Scan**' for leaf #5 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf #5 ($T_{\text{tar},5}$).

Configuration E: Transmittance Mode, Target

Step 62.

Carefully replace leaf #5 by leaf #6.

📌 NOTES

Etienne Laliberté 25 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces the inside of the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.

Configuration E: Transmittance Mode, Target

Step 63.

Collect a '**Target Scan**' for leaf #6 in this configuration and **save the file**.

📌 NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf #6 ($T_{\text{tar},6}$).

Configuration D: Transmittance Mode, Reference

Step 64.

Remove leaf #6 from the transmission sample port holder.

Configuration D: Transmittance Mode, Reference

Step 65.

Position leaf #1 over the **reflectance port** so that its **abaxial** (lower) side is now facing the inside of the sphere.

📌 NOTES

Etienne Laliberté 24 Apr 2018

Position the leaf to target the same area measured for the reflectance radiance, with the exception that its abaxial surface now faces the inside of the sphere. Position it so that the amount of leaf and vein material over the port is roughly proportional to the area of leaf and vein found throughout the leaf, while avoiding the large midrib vein. Position the leaf so that it is approximately halfway between the mid-rib vein and the leaf margin, and halfway between the tip and the base of the leaf lamina.


Configuration D: Transmittance Mode, Reference

Step 66.

Collect a '**Target Scan**' in this configuration and **save the file**.

Transmittance mode

Reference radiance (T_{ref})



NOTES

Etienne Laliberté 04 Jun 2018

This second reference radiance scan is only used to assess the stability of the system in transmittance mode.

Calculating Leaf Reflectance (Adaxial Surface)

Step 67.

The equation for **adaxial reflectance** of leaf i , $\rho_{\text{leaf},i}$ is

$$\rho_{\text{leaf},i} = [(R_{\text{tar},i} - R_{\text{str}}) \div (R_{\text{ref}} - R_{\text{str}})] \times \rho_{\text{ref}}$$

where

ρ_{ref} is the absolute reflectance of the calibrated Spectralon® reflectance standard.

Step 68.

The equation for **adaxial transmittance** of leaf i , $\tau_{\text{leaf},i}$ is

$$\tau_{\text{leaf},i} = (T_{\text{tar},i} \div T_{\text{ref}})$$

Warnings

The lamp of the integrating sphere can get **very hot** and should be handled from its slotted base to avoid burns.