

# General tube bacteria extraction (DNeasy)

Sarah Hessen-Schmidt

## Abstract

**Citation:** Sarah Hessen-Schmidt General tube bacteria extraction (DNeasy). **protocols.io**

dx.doi.org/10.17504/protocols.io.fbcbiw

**Published:** 13 Jul 2016

## Protocol

### Step 1.

Place Buffer ATL and Buffer AI at 56C for at least 5 minutes before starting. Make sure AW1 and AW2 have ethanol added

### Step 2.

Spin >500ul of bacteria liquid culture at 10000g for 10 minutes. Pour off supernatant

### Step 3.

Add 200ul buffer ATL to each sample and vortex for 15s

#### AMOUNT

200 µl Additional info:

#### REAGENTS

Buffer ATL (tissue lysis buffer) 19076 by [Qiagen](#)

#### DURATION

00:00:15

### Step 4.

Add 200ul 100% ethanol to each sample and vortex for 15s

#### AMOUNT

200 µl Additional info:

#### REAGENTS

✓ Ethanol by Contributed by users

#### DURATION

00:00:15

### Step 5.

Briefly spin to collect all liquid. Transfer all of sample including precipitate to spin column.

### Step 6.

Centrifuge at max speed (>20000g) for 2 minutes

 DURATION

00:02:00

### Step 7.

Empty collection tube

### Step 8.

Add 500ul buffer AW1 to spin column. Spin at max speed (>20000g) for 2 minutes

 AMOUNT

500 µl Additional info:

 REAGENTS

Buffer AW1 [19081](#) by [Qiagen](#)

 DURATION

00:02:00

### Step 9.

Add 500ul buffer AW2 to spin column. Spin at max speed (>20000g) for 2 minutes

 AMOUNT

500 µl Additional info:

 REAGENTS

Buffer AW2 [19072](#) by [Qiagen](#)

 DURATION

00:02:00

### Step 10.

Place spin column in a new collection tube and spin at max speed (>20000g) for 3 minutes

 DURATION

00:03:00

### Step 11.

Transfer spin column to a labeled 1.5ml Eppendorf lobind tube.

### Step 12.

Add 100ul of buffer TE directly to membrane/filter in spin column and incubate at room temp. for 1 min

 AMOUNT

100 µl Additional info:

 REAGENTS

✓ Buffer TE 1x by Contributed by users

🕒 DURATION

00:01:00

### Step 13.

Spin at 10000g for 2 minutes

🕒 DURATION

00:02:00