

# Minimal Inhibitory Concentration (MIC)

Brian Kenn Baltzar

## Abstract

This protocol is modified from Wiegand and colleagues protocol<sup>1</sup> to fit the conditions of our experiments.

<sup>1</sup>Wiegand, Irith, Kai Hilpert and Robert E W Hancock (2008). "Agar and broth dilution methods to determine the minimal inhibitory concentration (MIC) of antimicrobial substances". I: Nat Protoc 3.2, s. 163{75

**Citation:** Brian Kenn Baltzar Minimal Inhibitory Concentration (MIC). **protocols.io**

dx.doi.org/10.17504/protocols.io.gpwbvpe

**Published:** 14 Mar 2017

## Protocol

### Prepare medium

#### Step 1.

Müller-Hinton Agar and Broth

*Plate the agar on plates.*

### Prepare the bacterial isolates

#### Step 2.

Streak the bacteria onto agar plates without inhibitor.

#### Step 3.

Incubate plates for 18-24h at 37°C

### Prepare the bacterial isolates

#### Step 4.

For each isolate, select three to five morphologically similar colonies from the agar plates and transfer them to a glass tube with 5ml MH-Broth.

Incubate for 18-24h at 37°C

#### Step 5.

Measure OD450 on cultures from step 4, and make dilutions in x mL MHB to OD450 = 0,01

## Prepare the ELISA microplate

### Step 6.

Use the setup shown in the following figure.

	1	2	3	4	5	6	7	8	9	10	11	12	
A	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic A
B	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic B
C	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic C
D	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic D
E	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic E
F	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic F
G	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic G
H	X	X/2	X/4	X/8	X/16	X/32	X/64	X/128	X/256	X/512	Growth Control	Sterility Control	Antibiotic H

Growth control = Broth with bacteria inoculum but no antibiotic (growth)  
Sterility control = Only broth (no growth)

### Step 7.

Pipette 50µl of MHB into column 2-11

### Step 8.

Add 100µl MBH into column 12

### Step 9.

Add 100µl of the antibiotic into column 1

### Step 10.

Withdraw 50µl from each well in column 1, add this to the corresponding wells in column 2. Mix by pipetting up and down 4-6 times.

*(this makes column 2 a twofold dilution of column 1)*

### Step 11.

Repeat step 10 down to column 10 - discard the withdrawn solution from column 10.

## Inoculate the ELISA microplate

### Step 12.

Vortex the bacteria suspension from step 5.

Add 50µl of the suspension to each well in columns 1-11.

## Incubate the plates

### Step 13.

Cover the plates with film, make sure to create a tight seal to prevent any evaporation.

## Repeat

### Step 14.

To exclude any errors repeat every plate setup three times (triplicates).

## Incubate the plates

### Step 15.

Incubate the plate at 37°C for 16-24h.

## Results

### Step 16.

Measure the growth by determining OD<sub>450</sub>.

Compare the data from the three triplicates, and exclude any deviations