



Feb 27, 2019

Working

Neural progenitor banking

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x89frz6

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Celeste Karch

Washington University in St Louis

IPSC CORTICAL
DIFFERENTIATION
022017.pdf

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [IPSC CORTICAL DIFFERENTIATION](#) collection.

This method should be performed using sterile technique.

MATERIALS TEXT

Please refer to the attached full manuscript for required materials.

SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

- 1 Upon reaching at least 85% confluency, harvest neural progenitor cells as described in protocol below.
[Neural progenitor expansion protocol](#)
- 2 Perform a cell count in 3 ml of NIM using a hemacytometer.
- 3 Add equal volumes of NIM and 2x neural freezing medium to the NPC cell suspension for a final 1×10^6 cells/mL.
- 4 Gently mix solution and distribute 1 ml into sterile cryovials. Store cryovials in Styrofoam containers at -80 °C for 48:00:00 and then transfer to liquid nitrogen for long-term storage.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited