

# Gibson assembly

Ruihong Wang

## Abstract

Gibson is a PCR-based multi-fragment assembly technique. So whether PCR works has large impacts on final product.

**Citation:** Ruihong Wang Gibson assembly. **protocols.io**

dx.doi.org/10.17504/protocols.io.mzpc75n

**Published:** 15 Feb 2018


## Materials

🐛 Gibson Assembly Master Mix - 10 rxns [E2611S](#) by [New England Biolabs](#)

## Protocol

### Primer design

#### Step 1.



#### Step 2.

Construct a complete map of insert in the desired destination vectors in Snapgene or Ape.

### **Step 3.**

Identify the connection region between inserts and vector.

### **Step 4.**

Design primer for amplifying insert and vectors. Target melting temperature  $T_m$  set to 60 or higher.

### **Step 5.**

Order primers A, B, C and D.

*Notes: A has two parts. Starting from 5', the overhang is the reverse complement of Rev primer of backbone, and the primer region is a standard primer to amplify insert (Fwd primer of insert). B is same and has reverse complement of backbone Fwd primer followed by standard Rev primer of insert.*

So primer A= Backbone Rev primer (rev complement)+ Insert Fwd primer

Primer B= Backbone Fwd primer (rev complement)+ Insert Rev primer

Primer C is reverse complement of A.

Primer D is reverse complement of B.

### **Step 6.**

For more fragment to assemble, increase number of primers correspondingly.

## **PCR**

### **Step 7.**

Amplify insert with A and B primers. Amplify backbone with C and D primers. Use 25 or 50 uL reaction volume. Use DMSO or nested/ touchdown PCR if necessary.

### **Step 8.**

Add 0.25 uL DpnI and incubate at 37 C for 1hr (or longer depending on your paranoid) to remove template.

### **Step 9.**

Run a gel to verify the size of PCR product.

### **Step 10.**

Use PCR cleanup or Gel extraction to get clean product.

### **Step 11.**

Quantify the concentration by Nanodrop

## **Gibson Assembly**

### **Step 12.**

For 2 or 3 fragments assembly, use 0.02-0.5 pmols of DNA. For more fragments, use 0.2-1.0 pmols.

$\text{pmols} = (\text{weight in ng}) \times 1000 / (\text{base pairs} \times 650 \text{ daltons})$

**Step 13.**

Put insert and vectors into reaction and add DI water to 10 uL.

**Step 14.**

Add 10 uL 2X Gibson Master mix.

**Step 15.**

Incubate at 50 C for 15-60 minutes depending on the fragments number. 2-3 fragments use 15 mins. 4-6 use 60 mins.

**Step 16.**

Store in -20 C.

Transformation

**Step 17.**

Mix 1 uL of reaction and 10 uL of competent cells, do not vortex. Keep on ice for 30 mins.

**Step 18.**

Heat shock in 42 C water bath for 45 sec.

**Step 19.**

Keep on ice for 5 mins.

**Step 20.**

Add 300 uL SOC (pre-warm to Room temperature). Put in 250 rpm shaker in 37 C for 1hr.

**Step 21.**

Use 50-100 uL of recovered solution to plate out depending on copy number and plate size.

**Step 22.**

Harvest colony next day and pick single colony for liquid culture in corresponding antibiotics.

**Step 23.**

Sequencing and digestion to verify the assembly.