

# Lysate Preparation from Viruses in Cultured Cells Growing in a Monolayer

## Bioline

### Abstract

The steps for preparing the lysate are different depending on the starting material. Please ensure you follow the correct procedure for your starting material (see the section [Disrupting and Homogenizing Starting Materials](#) and the different lysate preparation protocols [here](#)). The subsequent steps detailed in [Genomic DNA removal and total RNA purification](#) are the same in all cases.

This protocol is for the isolation of integrated viral RNA from mammalian cells growing in culture.

**Citation:** Bioline Lysate Preparation from Viruses in Cultured Cells Growing in a Monolayer. **protocols.io**  
dx.doi.org/10.17504/protocols.io.f5dbq26

**Published:** 12 Dec 2016

## Guidelines

- A maximum of  $3 \times 10^6$  eukaryotic cells can be used as starting material. A hemocytometer can be used in conjunction with a microscope to count the number of cells. As a general guideline, a confluent 3.5 cm plate of HeLa cells will contain  $10^6$  cells.
- Cell pellets can be stored at  $-80^{\circ}\text{C}$  for later use or used directly in the procedure.  
Determine the number of cells present before freezing.
- Frozen pellets should be stored for no longer than 2 weeks to ensure that the integrity of the RNA is not compromised.
- Frozen cell pellets should not be thawed prior to beginning the protocol. Add the Lysis Buffer RX directly to the frozen cell pellet.

Additional reagents required: PBS (RNase-free).

Please review the Guidelines under [Genomic DNA removal and total RNA purification from all types of lysate](#) for other important details.

## Before start

- Ensure Lysis Buffer RX is prepared (see guidelines).
- Ensure that all solutions are at room temperature before use.
- Two types of spin columns are provided with this kit: the **Genomic DNA Removal Column** (blue ring) and the **RNA Column** (black ring). Ensure the correct column is used for each step of the procedure.
- All centrifugation steps are carried out in a benchtop microcentrifuge at 14,000 x g except where noted. Perform all centrifugation steps at room temperature.
- It is important to work quickly during this procedure.

## Materials

ISOLATE II Biofluids RNA Kit [BIO-52086](#) by [Bioline](#)

## Protocol

### Step 1.

Aspirate media and wash cell monolayer with an appropriate amount of PBS. Aspirate PBS.

#### NOTES

**Steve Hawkins** 08 Dec 2016

The steps for preparing the lysate are different depending on the starting material. Please ensure you follow the correct procedure for your starting material (see the different lysate protocols [here](#)). The subsequent steps (detailed [here](#)) for Genomic DNA Removal and Total RNA Purification are the same in all cases.

## Step 2.

Add 350 µL of Lysis Buffer RX directly to culture plate.

### 🔗 NOTES

**Steve Hawkins** 18 Oct 2016

Note: This volume of lysis buffer can be used for a culture plate of ~35 mm in diameter or up to one well of a 6-well culture dish.

## Step 3.

Lyse cells by gently tapping culture dish and swirling buffer around plate surface for 5 min.

### 🕒 DURATION

00:05:00

## Step 4.

Transfer lysate to a 1.5 mL RNase-free microcentrifuge tube (user supplied).

### 🔗 NOTES

**Steve Hawkins** 18 Oct 2016

Note: For input amounts greater than 10<sup>6</sup> cells, it is recommended that the lysate is passed through a nuclease-free 20 gauge (0.9 mm) syringe needle 5-10 times, in order to reduce the viscosity of the lysate prior to loading onto the column.

## Step 5.

Proceed to [Genomic DNA removal and total RNA purification from all types of lysate](#).

## Warnings

When working with chemicals, always wear a suitable lab coat, gloves and safety glasses.

Lysis Buffer RX contains guanidinium thiocyanate. This chemical is harmful in liquid form when in contact with skin or ingested. If the solution is allowed to dry, the powder is harmful if inhaled.

**CAUTION:** Do not add bleach directly to solutions or sample preparation waste containing guanidinium salts. Reactive compounds and toxic gases can form. In the case of spillage, clean the affected area with a suitable laboratory detergent and water.

For detailed information, please consult the material data safety sheet (MSDS) available on our website at [www.bioline.com](http://www.bioline.com).

Biofluids derived from all human and animal sources are considered potentially infectious. All necessary precautions recommended by the appropriate authorities in the country of use should be taken when working with biofluids.