



Aug 10,  
2019

# SPARC - Acute surgery and experimentation of the gastrointestinal tract and vagus nerve in the ferret

Charles C. Horn<sup>1</sup>, Derek M. Miller<sup>1</sup>, Stephanie Fulton<sup>1</sup>, Bill J. Yates<sup>1</sup>, Lee E. Fisher<sup>1</sup>, Ameya C. Nanivadekar<sup>1</sup>

<sup>1</sup>University of Pittsburgh

1 Works for me dx.doi.org/10.17504/protocols.io.6a7hahn

Charles C. Horn  
University of Pittsburgh

## ABSTRACT

This protocol is used for acute surgery and electrophysiological experimentation on the anesthetized ferret, specifically for manipulating the vagus nerve and gastrointestinal tract nerve stimulation and myoelectric recording.

Funding: This protocol was developed with funding from the NIH Common Fund's Stimulating Peripheral Activity to Relieve Conditions (SPARC) program (Award U18TR002205). To learn more about the SPARC program, visit <https://sparc.science>.

## EXTERNAL LINK

<https://www.biorxiv.org/content/10.1101/607242v2>

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

[Machine learning prediction of emesis and gastrointestinal state in ferrets Ameya C. Nanivadekar, Derek M. Miller, Stephanie Fulton, Liane Wong, John Ogren, Girish Chitnis, Bryan McLaughlin, Shuyan Zhai, Lee E. Fisher, Bill J. Yates, Charles C. Horn bioRxiv 607242; doi: <https://doi.org/10.1101/607242>](#)

## GUIDELINES

animal should be monitored for stability of respiration, heart rate, blood pressure, and temperature

## MATERIALS

NAME	CATALOG #	VENDOR
Vetbond	<a href="#">View</a>	3M corporation
30 ml barostat balloon	CT-BP-1017	MUI Scientific
3-point retractor	110-3445	George Tiemann

## Remove food

1 Remove food from the home cage 3 h before anesthesia 03:00:00

1.1 Measure body weight

## Preparation: Before first incision

2 Induce anesthesia using isoflurane in a chamber, then move to a facemask (1 to 2%).

- 3 Shave abdomen and ventral neck.
- 4 Place ferret supine on table. Connect EKG leads and place rectal probe for recording body temperature (keep at 36 to 40C).
- 5 Monitor heart rate and body temperature throughout procedure and testing.



CED Power 1401 running Spike 2 version 7.15  
electrophysiology recording  
CED CED Power 1401 running Spike 2 version 7.15

#### First incisions: Endotracheal tube and blood pressure catheter

- 6 Make a ~4cm incision in the ventral neck and dissection the trachea.
- 7 Insert an endotracheal tube and connect to pressure transducer to monitor respiration rate.



Small Animal Ventilator  
ventilator with expiration pressure transducer  
CWe SAR-830/AP

- 8 Make a ~3 cm incision on the medial surface of the right leg.
- 9 Place blood pressure catheter in left femoral artery.



Direct Blood Pressure System  
pressure transducer  
Kent Scientific DBP1000

#### Abdominal Surgery: Placement of electrodes and gastric tube/balloon

- 10 Make an incision in the midline abdominal muscle to expose the organs.
- 11 Open abdominal cavity using a 3-point retractor (Balfour pediatric retractor, George Tiemann), applied to the lateral edges of the abdominal incision and sternum. This retractor is elevated using a vertical post attached to a surgical board under the animal.

- 12 Make incision (~0.5 cm) on the left lateral edge of the gastric fundus and insert gastric catheter with the tip resting in the antrum; in parallel with the gastric catheter is a 30 ml barostat balloon. Place a purse-string suture around the incision and tied. And, additional knot is made around the catheter and the free ends of this knot and purse-string suture are tied together. Finally surgical glue (Vetbond) is applied to the incision area.
- 13 In some experiments, a cuff electrode is placed on either the ventral or dorsal abdominal vagus nerve trunk. Suture a ground electrode lead to connective tissue on the greater curvature of the stomach.
- 13.1 This is accomplished by caudal retraction of the stomach and retraction of the liver to the right side of the animal using saline moistened gauze. In this position, the esophagus and vagi become accessible.
- 14 Place planar electrodes on surface of stomach and duodenum using 8-0 silk suture. Suture a ground electrode lead to connective tissue on the greater curvature of the stomach
- 15 Test electrode impedances.



Grapevine Neural Interface Processor and a Nano2 recording headstage electrophysiological recording and stimulation  
Ripple LLC Grapevine Neural Interface Processor and a Nano2 r

#### Testing

- 16 Following this initial setup various tests can be performed, including testing the effects of stomach distension, nerve electrical stimulation, and emetine infusion on gastrointestinal myoelectric responses. Emetine is a chemical emetic agent that produces retching and emesis when infused into the stomach and small intestine. Because emesis can produce strong abdominal contractions it is recommended that the abdominal incision is securely sutured.

#### Gastric distension trials

- 17 Measure impedances of electrodes. Record 5 minutes baseline before and after each distension trial.



Grapevine Neural Interface Processor and a Nano2 recording headstage electrophysiological recording and stimulation  
Ripple LLC Grapevine Neural Interface Processor and a Nano2 r

- 18 Fill balloon using a syringe pump with 5 ml, 10 ml, or 20 ml saline. Infusion rate is 10ml/min. Saline is held in balloon for 5 minutes before deflation.

## Emetine Trial

- 19 Record 20 min of baseline gastrointestinal myoelectric activity.



Grapevine Neural Interface Processor and a Nano2 recording headstage electrophysiological recording and stimulation

Ripple LLC Grapevine Neural Interface Processor and a Nano2 r

- 20 Use a syringe pump to infuse 5 mg/kg emetine into the stomach. Record for 1 hour.

## Cuff Stim Trials

- 21 Apply desired stimulation parameters (pulse width, current amplitude, frequency) to abdominal vagus while recording. Record baseline data before and after trials and measure impedances at the end.



Grapevine Neural Interface Processor and a Nano2 recording headstage electrophysiological recording and stimulation

Ripple LLC Grapevine Neural Interface Processor and a Nano2 r

- 22 Open abdomen and remove electrodes, balloon, and catheter.

- 23 Euthanize animal with intracardiac euthasol injection.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited