

# SypHy imaging and dextran uptake assays from Yao CK et al. (2017)

Chi-Kuang Yao, Yu-Tzu Liu, I-Chi Lee, You-Tung Wang, Ping-Yen Wu

## Abstract

This protocol is from 'Flower  $\text{Ca}^{2+}$  channel in CME and ADBE' of Yao CK et al.

Please see the manuscript for the full method details.

**Citation:** Chi-Kuang Yao, Yu-Tzu Liu, I-Chi Lee, You-Tung Wang, Ping-Yen Wu SypHy imaging and dextran uptake assays from Yao CK et al. (2017). **protocols.io**

[dx.doi.org/10.17504/protocols.io.hhbb32n](https://doi.org/10.17504/protocols.io.hhbb32n)

**Published:** 04 Apr 2017

## Guidelines

Experiments were performed at either room or physiological temperatures controlled by Warner temperature controller (TC-344B).

## Before start

You'll need:

**chamber (Warner instruments RC-25F) with two parallel platinum wires separated by 5 mm**

### imaging buffer:

- 136 mM NaCl
- 2.5 mM KCl
- 2 mM  $\text{CaCl}_2$
- 1.3 mM  $\text{MgCl}_2$
- 10 mM glucose
- 10 mM HEPES [pH 7.4]
- 10 mM CNQX
- 50 mM AP-5

### acidic buffer:

- 136 mM NaCl
- 2.5 mM KCl
- 2 mM  $\text{CaCl}_2$
- 1.3 mM  $\text{MgCl}_2$
- 10 mM glucose
- 10 mM 2-[N-morpholino] ethane sulphonic acid [pH 5.5]
- 10 mM CNQX
- 50 mM AP-5

Acidic buffer was prepared by replacing HEPES in the imaging buffer with 2-[N-morpholino] ethane sulphonic acid [[47](#),[48](#)].

### imaging solution:

- 144 mM NaCl
- 2.5 mM KCl
- 2.5 mM  $\text{CaCl}_2$
- 2.5 mM  $\text{MgCl}_2$
- 10 mM HEPES (pH 7.5)
- 10 mM CNQX
- 50 mM AP-5)

## Materials

- ✓ imaging buffer by Contributed by users
- ✓ acidic buffer by Contributed by users
- ✓ imaging solution by Contributed by users
- ✓ pSpCas9(BB) by Contributed by users
- ✓ pCMV-SyphyA4 #24478 by [addgene](#)
- ✓ 50  $\mu\text{M}$  40 kDa TMR-dextran (Invitrogen) by Contributed by users

## Protocol

### SypHy imaging

#### Step 1.

Transfect DIV7 cultured rat hippocampal neurons cultured on 12-well plate with pSpCas9(BB) and pCMV-SyphyA4 (addgene#24478) plasmids [5] by  $\text{Ca}^{2+}$  phosphate method.

### SypHy imaging

#### Step 2.

Bath DIV13-15 neurons in the imaging buffer in a chamber (Warner instruments RC-25F) with two parallel platinum wires separated by 5 mm.

### SypHy imaging

#### Step 3.

Elicit SV exocytosis with a train of 200 action potentials delivered with 20-Hz electric field stimulation (50 mA, 1-ms pulse width).

#### 📌 NOTES

**Chi-Kuang Yao** 04 Apr 2017

See "before start" for the imaging buffer.

### SypHy imaging

#### Step 4.

Capture single images every 1 s using MetaMorph software and ANDOR iXon 897 camera.

### SV re-acidification

#### Step 5.

For SV re-acidification experiments, perfuse the imaging chamber with the imaging buffer.

### SV re-acidification

#### Step 6.

Perfuse with an acidic buffer (see 'Before start' for exact composition).

### SV re-acidification

#### Step 7.

Perfuse the imaging buffer to allow surface SypHy to be fluorescent.

#### 📌 NOTES

**Chi-Kuang Yao** 29 Mar 2017

Maintain experimental temperatures at physiological temperatures.

## Counting

### Step 8.

Count the final SypHy fluorescence intensity in the presynaptic terminals by subtracting the background fluorescence intensity on the surrounding coverslip from the SypHy fluorescence intensity within presynaptic terminals.

#### 📌 NOTES

**Chi-Kuang Yao** 29 Mar 2017

Each data value was obtained from a single terminal.

## 40 kDa TMR-dextran uptake assays

### Step 9.

For 40 kDa TMR-dextran uptake assays, stimulate DIV13-15 neurons transfected with pSpCas9(BB) plasmids by a train of 1600 action potentials delivered with a 80 Hz electric field stimulation (50 mA, 1-ms pulse width) in the imaging solution in the presence of 50  $\mu$ M 40 kDa TMR-dextran (Invitrogen).

## 40 kDa TMR-dextran uptake assays

### Step 10.

Perfuse neurons with the same buffer for 5 min to remove excess dextran dye.

#### 🕒 DURATION

00:05:00

#### 📌 NOTES

**Chi-Kuang Yao** 29 Mar 2017

Experiments were performed at room temperature. Imaging was achieved through MetaMorph software and ANDOR iXon 897 camera. Image processing was achieved using Image J and LSM Zen.