



May 17,  
2019

Working

## Yale - Blood and Urine Inorganic Phosphorous [↗](#)

John Stack<sup>1</sup>, Gary Cline<sup>1</sup>

<sup>1</sup>Yale University

[dx.doi.org/10.17504/protocols.io.y2gfybw](https://doi.org/10.17504/protocols.io.y2gfybw)

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)

Lili Liang

### ABSTRACT

#### Summary:

Procedure used to determine the concentration of inorganic phosphorus in blood, serum, and plasma. Phosphorus (Inorganic) is measured as the complex with molybdate and monitored at 340nm.

### EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=217&docType=Protocol>

### MATERIALS

NAME	CATALOG #	VENDOR
Phosphorus Reagent	R85562	Prolabs(cliniqa)
Multi Analyte Calibrator	R60010	Prolabs(cliniqa)
Assayed Control Serum 1	R83082	Prolabs(cliniqa)
Assayed Control Serum 2	R83083	Prolabs(cliniqa)

### MATERIALS TEXT

#### Reagent Preparation:

**Phosphorus Reagent:** As supplied by vendor

**Multi Analyte Calibrator:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

**Assayed Control Serum 1:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

**Assayed Control Serum 2:** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

### BEFORE STARTING

*Analysis by automated system Cobas Mira Plus*

- 1 Calibrate Cobas for Phosphorus analysis by running a multi-analyte serum calibrator and two control serum.

- 2 Sample handling as performed by Cobas Mira Plus
- Pipette 3  $\mu\text{L}$  of sample into cuvette.
  - Add 280  $\mu\text{L}$  of Phosphorus Liquid Reagent.
  - Mixture is incubated at 37°C for 10 minutes
  - Absorbance is measured at 340nm.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited