

Protocol F: Preparation of recovery buffer

David Booth

Abstract

Citation: David Booth Protocol F: Preparation of recovery buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.h82b9ye

Published: 22 Jun 2018

Protocol

Recovery buffer recipe

Step 1.

10 mM HEPES-KOH, pH 7.5, 0.9 M D-sorbitol, 8% (w/v) PEG-8000

Sterile filter buffer through 0.22 µm filter

Step 2.

Store at 4°C

Step 3.