

Working

## Preparation of nasopharyngeal samples for immunofluorescence assay 👄

PLOS One

Maria Lucia Silva<sup>1</sup>, Flavia Escremim de Paula<sup>1</sup>, Talita Gagliardi<sup>2</sup>, Eurico Arruda<sup>1</sup>

<sup>1</sup>School of Medicine - University of Sao Paulo, <sup>2</sup>University of Maryland, College Park dx.doi.org/10.17504/protocols.io.w8rfhv6

Talita B. Gagliardi
University of Maryland, College Park

**EXTERNAL LINK** 

## https://doi.org/10.1371/journal.pone.0217744

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Matsuno AK, Gagliardi TB, Paula FE, Luna LKS, Jesus BLS, Stein RT, Aragon DC, Carlotti APCP, Arruda E (2019) Human coronavirus alone or in co-infection with rhinovirus C is a risk factor for severe respiratory disease and admission to the pediatric intensive care unit: A one-year study in Southeast Brazil. PLoS ONE 14(6): e0217744. doi: 10.1371/journal.pone.0217744

- 1 Add the nasopharyngeal sample into a conical tube.
- 2 Only in samples with much mucus, drop N-L-Acetyl-Cysteine (maximum 1ml) into the tube and homogenize by vortex.
- 3 Centrifuge the tube for 10 minutes at 4oC and 300 x g.
- Discard the supernatant.
- 5 Homogenize the pellet of cells with phosphate-buffered saline solution 1x (PBS 1x).
- 6 Drop 5ul of diluted suspension of cells into a glass slide.
- 7 Let the slide inside the hood or in front of a hair dryer until the sample dries.
- 8 Immerge the glass slide into pure and cold acetone (fixing solution).

- 9 After 5 minutes, wash the slide one in PBS 1x.
- 10 Store the glass slides at -20oC.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited