

30.0% Acrylamide stock solution (30.0:0.8)

David Dunigan and Irina Agarkova

Abstract

For use in [Polyacrylamide Gel System For Electrophoresis Of Proteins](#).

Citation: David Dunigan and Irina Agarkova 30.0% Acrylamide stock solution (30.0:0.8). **protocols.io**
dx.doi.org/10.17504/protocols.io.eq2bdye

Published: 13 Jun 2016

Protocol

Step 1.

Add 60.0 gm of acrylamide and 1.6 gm of N, N'-methylenebisacrylamide (bisacrylamide) to 140 mL of d-H₂O.

Step 2.

Stir until dissolved.

Step 3.

Adjust the volume to 200 mL.

Step 4.

Store protected from the light.