

Bacterial transformation version 2

Stephen Floor

Abstract

Standard protocol to transform bacteria with a plasmid using chemically competent E. coli and antibiotic resistance.

Citation: Stephen Floor Bacterial transformation. **protocols.io**

dx.doi.org/10.17504/protocols.io.mnwc5fe

Published: 16 Jan 2018

Materials

- ✓ MACH1 or DH5a or TOP10 or NEB stable or etc for DNA purification/miniprep by Contributed by users
- ✓ BL21(DE3) or BL21-Star(DE3) or Rosetta2(DE3) or etc for protein purification by Contributed by users

Protocol

Step 1.

Count the number of transformations to perform and remove the appropriate number of aliquots from the -80 to an ice bucket. Our aliquots are 100 ul each.

🔗 NOTES

Stephen Floor 16 Jan 2018

If transforming a plasmid with repeats (such as a lentiviral plasmid), be sure to use NEB Stable or Stbl3 cells or similar.

Step 2.

Thaw cells on ice.

🔗 NOTES

Stephen Floor 16 Jan 2018

All steps are in a microcentrifuge tube.

Step 3.

Label microcentrifuge tubes and add 20 ul cells to each tube.

Step 4.

Add 10-500 ng DNA (typically just use 1 uL) to each tube of cells.

Step 5.

Incubate for 25 minutes on ice. Anywhere from 5 to 25 minutes usually works.

Step 6.

Heat shock at 42 degrees for one minute in water bath.

 [TEMPERATURE](#)

42 °C Additional info:

Step 7.

Recover for two minutes on ice.

Step 8.

Add 180 µL LB or SOC media, mix.

 [AMOUNT](#)

180 µl Additional info: LB or SOC media

Step 9.

Grow at 37 degrees for one hour with shaking. Warm the plates during this step.

 [TEMPERATURE](#)

37 °C Additional info:

 [NOTES](#)

Stephen Floor 16 Jan 2018

If a faster transformation is desired the growth step can be skipped entirely for ampicillin resistant plasmids, or shortened to ~15 minutes for kanamycin resistant plasmids. Growth may increase the number of colonies though.

Step 10.

Plate 75 µL of the transformation on plate with appropriate antibiotic.

 [ANNOTATIONS](#)

Stephen Floor 16 Jan 2018

If the transformation is low efficiency (cloning), plate as much as 200 ul.

Step 11.

Leave at 37 degrees overnight or room temperature for two days.