

# 1% SDS in DNA Buffer

Andrew Baker

## Abstract

This buffer is used to stabilize samples for archiving and subsequent genomic DNA extraction.

SDS is dissolved to a concentration of 1% (w/v) in DNAB (DNA buffer: 0.4 M NaCl + 0.05 M EDTA in MilliQ water). The buffer may need to be warmed for SDS to completely dissolve. Tissue samples are added to small aliquots of 1% SDS in DNAB and heated to 65°C for 60-90 minutes. They are then stable at room temperature and ready for extraction of genomic DNA.

**Citation:** Andrew Baker 1% SDS in DNA Buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.dys7wd

**Published:** 04 Mar 2016

## Guidelines

SDS is dissolved to a concentration of 1% (w/v) in DNAB (DNA buffer: 0.4 M NaCl + 0.05 M EDTA in MilliQ water). The buffer may need to be warmed for SDS to completely dissolve. Tissue samples are added to small aliquots of 1% SDS in DNAB and heated to 65°C for 60-90 minutes. They are then stable at room temperature and ready for extraction of genomic DNA.

## Protocol

### Step 1.

Prepare stock solution of 4 M Sodium chloride in MilliQ water



#### REAGENTS



Sodium chloride [View](#) by [P212121](#)

### Step 2.

Prepare stock solution of 0.5 M EDTA in MilliQ water



#### REAGENTS



Ethylenediaminetetraacetic acid by Contributed by users

### Step 3.

Mix 50 mL 4 M NaCl and 50 mL 0.5 M EDTA

### Step 4.

Make up to a final volume of 500 mL with MilliQ water



#### REAGENTS



MilliQ water by Contributed by users

### Step 5.

Dissolve SDS in DNA Buffer to a final concentration of 1% (w/v).  
e.g., 5 g SDS in 500 mL of DNA Buffer.



## REAGENTS



Sodium Dodecyl Sulfate [View](#) by [P212121](#)