

# pH tolerance assay for *Vibrio natriegens*

Daniel Marchal

## Abstract

**Citation:** Daniel Marchal pH tolerance assay for *Vibrio natriegens*. **protocols.io**

[dx.doi.org/10.17504/protocols.io.prjdm4n](https://dx.doi.org/10.17504/protocols.io.prjdm4n)

**Published:** 25 Apr 2018

## Before start

Media preparation:

1M  $\text{KH}_2\text{PO}_4$  100ml → 13,61g

1M  $\text{K}_2\text{HPO}_4$  100ml → 17,42g

1M Na-Acetate 100ml → 13,61g Trihydrate

1M Acetic acid 100ml → 60,05g (liquid)

1M Glycine 100ml → 7,51g

Autoclave  $\text{KH}_2\text{PO}_4$ ,  $\text{K}_2\text{HPO}_4$  and Na-Acetate and sterilfiltrate Glycine, Acetic acid should be sterile yet

LB2,5:

25g/l LB broth (Miller)

15g/l NaCl

BHlv2:

37g/l Brain heart infusion broth

204 mM NaCl

4.2 mM KCl

23.14 mM MgCl<sub>2</sub>

## Materials

● Sodium acetate [View](#) by [P212121](#)

● Potassium phosphate (dibasic) [View](#) by [P212121](#)

✓ Lysogeny broth by Contributed by users

✓ Hydrochloric Acid by Contributed by users

✓ brain Heart Infusion Broth Oxoid CM1135-UK by Contributed by users

Sodium hydroxide [SB0617.SIZE.500G](#) by [Bio Basic Inc.](#)

✓ Monopotassium phosphate by Contributed by users

Acetic acid, glacial 1.01830.2500 by [Merck Millipore](#)

Glycine [50046](#) by [Sigma](#)

## Protocol

### Step 1.

Inoculate precultur of BHIv2 with *V. natriegens* and incubate oN at 30°C

### Step 2.

Prepare buffer solutions for each pH value in the table.

All stock solutions should be at the same molarity (volumens refer to an buffer volume of 50ml):

pH	Acid	V( Acid-Stock)	Base	V(Base-Stock)
2.4	HCl	32,40	Glycine	50,00
3	HCl	11,40	Glycine	50,00
3.5	Acetic acid	2,60	Sodium acetate	47,40
4	Acetic acid	7,40	Sodium acetate	42,60
4.5	Acetic acid	17,73	Sodium acetate	32,27
5	Acetic acid	31,74	Sodium acetate	18,26
5.5	Acetic acid	42,30	Sodium acetate	7,70
6	KH <sub>2</sub> PO <sub>4</sub>	2,90	K <sub>2</sub> HPO <sub>4</sub>	47,10
6.5	KH <sub>2</sub> PO <sub>4</sub>	8,16	K <sub>2</sub> HPO <sub>4</sub>	41,84
7	KH <sub>2</sub> PO <sub>4</sub>	19,07	K <sub>2</sub> HPO <sub>4</sub>	30,93
7.5	KH <sub>2</sub> PO <sub>4</sub>	33,05	K <sub>2</sub> HPO <sub>4</sub>	16,95
8	KH <sub>2</sub> PO <sub>4</sub>	43,02	K <sub>2</sub> HPO <sub>4</sub>	6,98
8.5	Glycine	50,00	NaOH	4,00
9	Glycine	50,00	NaOH	8,80

9.5	Glycine	50,00	NaOH	22,40
-----	---------	-------	------	-------

### **Step 3.**

Measure pH and adjust with HCl/NaOH

### **Step 4.**

Prepare 5ml flasks with 4ml BHI + 1ml buffer (see table) for pH2 – pH9

### **Step 5.**

Measure pH and adjust with HCl/NaOH

### **Step 6.**

Sterilfiltrate into a new tube (because the medium is contaminated by the pH meter)

### **Step 7.**

Inoculate flasks with 10µl preculture

### **Step 8.**

Mix well and incubate at 37°C shaking

### **Step 9.**

Take samples after 15min and 6h and make 1:10, 1:100 and 1:1000 dilutions

### **Step 10.**

Plate out 50µl of each dilution on LB2,5 medium

### **Step 11.**

Incubate plates at 37°C oN and determine CFUs