



May 10,
2019

Working

U Mass - Uric Acid [↗](#)

Jason Kim¹

¹University of Massachusetts

dx.doi.org/10.17504/protocols.io.x5sfq6e

Mouse Metabolic Phenotyping Centers

Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

This experiment involves a spectrophotometric measurement using Roche Cobas Clinical Chemistry Analyzer. Serum uric acid levels are affected by alterations in systemic protein and nitrogen metabolism. Serum uric acid levels are altered in kidney failure and renal complications of diabetes.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=188&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
Calibrator f.a.s.	10759350 360	Roche		
Precinorm U plus	12149435 160	Roche		
Precipath U plus	12149443 160	Roche		
NaCl Diluent 9 %	04774230 190	Roche		
Chimneys	11930630 001	Roche		
Cleaner	04774248 190	Roche		
Micro Sample cups	11406680 001	Roche		
NERL High Quality Water	9805	Fisher Scientific		
Uric Acid	04657608190	Roche		

MATERIALS TEXT

Note:

Roche, [RRID:SCR_001326](#)

Fisher Scientific, [RRID:SCR_008452](#)

BEFORE STARTING

Notes:

✓ Try to use freshly prepared serum and plasma samples for this assay.

✓ No dilution or treatment of the sample is required, but plasma samples should be centrifuged to remove any fibrin/fibrinogen clumps.

✓ Samples should be stored at 2-8°C for 24 hours prior to analysis. For longer periods, store samples at -70°C, and avoid repeated freeze/thaw cycles.

✓ A 50 µl dead volume is required in addition to sample volume for multi-protein analysis (typically 1-5 µl).

- 1 Perform daily quality control assessment of instrumentation before analysis.
- 2 Load each sample into a specialized micro-sample cup for the clinical chemistry analyzer.
- 3 Select Uric acid test on display and run the analysis.
- 4 Collect and analyze the data.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited