

# Protocol B: Preparation of ATP for transfection

David Booth

## Abstract

**Citation:** David Booth Protocol B: Preparation of ATP for transfection. **protocols.io**

dx.doi.org/10.17504/protocols.io.h8ub9ww

**Published:** 22 Jun 2018

## Protocol

Dissolve solid disodium adenosine triphosphate (Sigma A7699) in a minimal volume of sterile water

### Step 1.

Titrate ATP solution with 1 M NaOH to pH 7.5, monitoring the pH with pH strips

### Step 2.

Sterile filter ATP solution through a 0.22  $\mu\text{m}$  filter

### Step 3.

Quantify ATP using UV-Vis spectrophotometry

### Step 4.

- Prepare a small aliquot of a 100 fold dilution of ATP in water
- Measure the optical absorbance at 259 nm with a UV-Vis spectrophotometer
- Determine the concentration of ATP with an extinction coefficient of  $15,400 \text{ M}^{-1}\text{cm}^{-1}$ .

Dilute ATP to a final concentration of 250 mM

### Step 5.

- Adjust volume from step 4 with water that has been sterily filtered through a 0.22  $\mu\text{m}$  filter

Store at  $-20^{\circ}\text{C}$

### Step 6.