

# Staining Unfixed Labyrinthulomycetes with Nile Red

This protocol modified from: Pandey A, Bhathena Z. 2014. Prevalence of PUFA Rich Thraustochytrids sps. along the Coast of Mumbai for Produ

## Abstract

This protocol modified from:

Pandey A, Bhathena Z. 2014. Prevalence of PUFA Rich Thraustochytrids sps. along the Coast of Mumbai for Production of Bio Oil. Journal of Food and Nutrition Research 2(12): 993-999.

**Citation:** This protocol modified from: Pandey A, Bhathena Z. 2014. Prevalence of PUFA Rich Thraustochytrids sps. along the Coast of Mumbai for Produ Staining Unfixed Labyrinthulomycetes with Nile Red. **protocols.io** dx.doi.org/10.17504/protocols.io.hghb3t6

**Published:** 30 Mar 2017

## Protocol

### Step 1.

Primary Nile Red (NR) stock: 1mg/mL in DMSO.

Working NR stock: 100 µg/mL in DMSO.

\*Keep NR in the dark.\*

\*NR is not very soluble in water\*

We have used this protocol successfully for *Schizochytrium* 28209 grown in 1/2 790 and *Artemia* soup.

We've even seen stained zoospores still swimming!

### 📌 NOTES

**Laura Halligan** 30 Mar 2017

**ATCC Medium: 790 By+ Medium**

Yeast Extract.....1.0 g

Peptone.....1.0 g

D+-Glucose.....5.0 g

Seawater.....1000 ml

<https://www.atcc.org/~media/920FDAC93FF84B79851C29FBB8049862.ashx>

**'Artemia Soup' medium**

1.25 g ground, freeze-dried brine shrimp

autoclaved in 1 liter artificial seawater

### Step 2.

Add 10 µl of NR working stock to 100 µl of cell culture.

Final staining condition: 10 µg/mL in 10% DMSO.

### **Step 3.**

Cover test tube (wrap in aluminum foil) to prevent NR light exposure and then vortex for 1 minute.

### **Step 4.**

Let stained sample incubate in the dark at room temperature for 5 minutes.

### **Step 5.**

After 5 minutes, mount cells on slide and place cover slip.

### **Step 6.**

The photo shown was taken with a Lietz I3 filter cube (excitation BP 450-490 nm; dichroic 510 nm; emission LP 515).

### **Step 7.**