

T-maze phototaxis with LEDs

Björn Brembs

Abstract

We use this protocol to test if flies have preference for darkness or light (phototactic or non-phototactic response).

Citation: Björn Brembs T-maze phototaxis with LEDs. **protocols.io**

dx.doi.org/10.17504/protocols.io.c76zrd

Published: 28 Jul 2015

Guidelines

We are using 3 day-old-flies (half of the experimental flies have their wings clipped, and half don't).

Protocol

Step 1.

Clip wings of half of the flies.

NOTES

Björn Brembs 20 Jul 2015

We are using 3 day-old-flies (half of the experimental flies have their wings clipped, and half don't).

Step 2.

Let flies recover for 1 day (at 25°C, faster recovery).

DURATION

24:00:00

Step 3.

Put the flies in the T-maze (entrance tube) and let them adapt themselves to the machine and the darkness for 10 minutes.

DURATION

00:10:00

Step 4.

Move the flies into the elevator (the mobile part of the T-maze), and let them stay there for around 1 minute.

DURATION

00:01:00

Step 5.

Move the flies to the bottom of T-maze (taking down the elevator).

Now they will have to make a decision between two different tubes: one is completely dark, and the other one has a LED in the bottom/is entirely transparent.

We are going to use different LEDs:

White LEDs
Yellow LEDs
Green LEDs

Step 6.

Wait for 30 seconds for the flies to make their decision, and lock them in the tubes that they have chosen (taking up the elevator).

 DURATION

00:00:30

Step 7.

Anesthetize the flies with CO₂ and count them under the magnifying glass.