

Large (Klenow) Fragment Blunting (M0210)

New England Biolabs

Abstract


Protocol for blunting ends by 3' overhang removal and fill-in of 3' recessed (5' overhang) ends using DNA Polymerase I, Large (Klenow) Fragment (M0210)

Citation: New England Biolabs Large (Klenow) Fragment Blunting (M0210). [protocols.io](https://www.protocols.io)

dx.doi.org/10.17504/protocols.io.cpgvjv

Published: 23 Jan 2015

Materials

 DNA Polymerase I Klenow Fragment - 200 units [M0210S](#) by [New England Biolabs](#)

Protocol

Step 1.

DNA should be dissolved in 1x NEBuffer or T4 DNA Ligase Reaction buffer supplemented with 33 μ M each dNTP.

ANNOTATIONS


Alexander Chamessian 18 Dec 2015

I went right into a BstAPI digest that was done in Cutsmart. NEB says that DNA Poly Klenow has 100% activity in Cutsmart buffer.

Step 2.

Add 1 unit of Klenow per microgram DNA.

REAGENTS

 DNA Polymerase I Klenow Fragment - 200 units [M0210S](#) by [New England Biolabs](#)

Step 3.

Incubate for 15 minutes at 25°C.

DURATION

00:15:00

Step 4.

Stop reaction by adding EDTA to a final concentration of 10 mM.

ANNOTATIONS

Alexander Chamessian 18 Dec 2015

Rather than EDTA, I just went for the heat inactivation, using 80C instead.

Step 5.

Heat for 20 minutes at 75°C.

DURATION

00:20:00