OPEN ACCESS



On-Column DNase I Treatment for ISOLATE II Biofluids RNA Kit

Bioline

Abstract

The ISOLATE II Biofluids RNA Kit isolates total RNA with minimal amounts of genomic DNA contamination using the supplied Genomic DNA Removal Column. However, additional DNase I treatment may be required in certain cases e.g. the amount of genomic DNA in the sample exceeds the capacity of the Genomic DNA Removal Column, or performing a highly sensitive application.

An alternative additional DNase I treatment protocol is: <u>DNase I Treatment of Purified RNA in Solution</u>.

This optional protocol can be used for additional removal of residual DNA that may affect sensitive downstream applications.

Citation: Bioline On-Column DNase I Treatment for ISOLATE II Biofluids RNA Kit. protocols.io

dx.doi.org/10.17504/protocols.io.f5ibg4e

Published: 12 Dec 2016

Guidelines

Please review the Guidelines under <u>Genomic DNA removal and total RNA purification from all types of</u> lysate for important details.

Materials

ISOLATE II Biofluids RNA Kit BIO-52086 by Bioline

Protocol

Step 1.

For each on-column digest to be performed, prepare a DNase I - buffer mix by adding 1 5μ L of the supplied DNase I Solution to 100 μ L of DNase I Reaction Buffer DRB.

Step 2.

Mix gently by inverting the tube a few times. Do not vortex.

Step 3.

Perform the appropriate RNA isolation procedure for your starting material up to and including the 'Binding RNA to Column' section (up to step 12 of <u>Genomic DNA removal and total RNA purification from all types of lysate</u>).

Step 4.

Apply 400 µL of Wash Buffer W1 to the column and centrifuge for 2 min at 14,000 x g.

O DURATION

00:02:00

Step 5.

Discard the flow-through. Reassemble the spin column with its Collection Tube.

P NOTES

Steve Hawkins 18 Oct 2016

Note: Ensure the entire wash buffer volume has passed through into the Collection Tube by inspecting the column. If the entire wash volume has not passed, spin for an additional minute at $14,000 \times g$.

Step 6.

Apply 115 µL of the DNase I - buffer mix to the column and centrifuge for 1 min at 14,000 x g.

O DURATION

00:01:00

NOTES

Steve Hawkins 18 Oct 2016

Note: Ensure the entire volume of DNase I - buffer mix passes through the column. If necessary, spin for an additional minute at $14,000 \times g$.

Step 7.

Pipette the flow-through present in the Collection Tube back on to the top of the column.

NOTES

Steve Hawkins 18 Oct 2016

Note: This step must be performed to ensure maximum DNase I activity and to obtain maximum yields of RNA. This is particularly important for the isolation of small RNA species.

Step 8.

Incubate at room temperature (18-25°C) for 15 min.

© DURATION 00:15:00

Step 9.

Without any further centrifugation, proceed directly to the second wash step in the RNA Column Wash section (step 15 here). Apply the wash buffer directly to the column containing the DNase I - buffer mix.