

# Dot Blot Preparation

David Dunigan and Irina Agarkova

## Abstract

**Citation:** David Dunigan and Irina Agarkova Dot Blot Preparation. **protocols.io**

[dx.doi.org/10.17504/protocols.io.etibeke](https://dx.doi.org/10.17504/protocols.io.etibeke)


**Published:** 13 Jun 2016

## Guidelines

### Materials:

- 1.) host and viral DNAs
- 2.) 1X TE buffer (10 mM Tris, pH 8.0, 1 mM EDTA)
- 3.) 1 M NaOH
- 4.) 1 M Acetic acid
- 5.) nylon membrane, cut to 12.3 X16 cm (2 sheets)
- 6.) Whatman 3MM filter paper, cut 12.3X16 cm (2 sheets)
- 7.) SSC solutions: 1X, 5X, 20X  
20X SSC (1000 mL):
  - 175.3 g NaCl
  - 88.2 g sodium citrate
  - adjust pH to 7.0 with NaOH
- 8.) 96-well assay plates
- 9.) hybridot dot blot apparatus and pump

Note: Load the membranes so that when scanned, the second blot of the set may be placed under the first, i.e., load as the following diagram shows:



## Protocol

### DNA dilutions

#### Step 1.

Prepare DNAs in 1.2 µg/48.6 µl concentrations.

### DNA dilutions

#### Step 2.

Distribute 97.2 µl of each to 96-well plates in the same pattern as the 1.0 µg samples will be loaded onto the membrane.

### DNA dilutions

#### Step 3.

Make 3 serial 1/2 dilutions with each DNA sample (48.6 µl into 48.6 µl 1XTE).

### DNA dilutions

#### Step 4.

Discard 48.6 µl from the final dilution.

### DNA dilutions

#### Step 5.

Add 5.4 µl of 1 M NaOH to each sample.

### DNA dilutions

#### Step 6.

Cover plates with parafilm.

### DNA dilutions

#### Step 7.

Incubate at 37°C, 30 min.

 DURATION

00:30:00

### DNA dilutions

#### Step 8.

Add 6.0 µl of 1 M acetic acid to each sample.

### DNA dilutions

#### Step 9.

Chill on ice.

### Membrane preparation

#### Step 10.

Wet the membrane in 1X SSC ≥60 min.

 DURATION

01:00:00

### Membrane preparation

#### Step 11.

On the hybridot manifold, place a piece of Whatman 3MM filter paper and wet with 20X SSC.

### Membrane preparation

#### Step 12.

Roll air pockets out.

### Membrane preparation

#### Step 13.

Mount the membrane on top of the 3MM paper (centered).

#### Membrane preparation

##### **Step 14.**

Roll air pockets out.

#### Membrane preparation

##### **Step 15.**

Wash the membrane with 5X SSC (ca. 25 ml) with vacuum filtration.

#### Membrane preparation

##### **Step 16.**

Screw the top of the hybridot manifold in place, then shut off the vacuum pump.

#### Membrane preparation

##### **Step 17.**

Add 200 µl of 5X SSC to each well.

#### Membrane preparation

##### **Step 18.**

Add DNAs (50 µl each) to the wells.

#### Membrane preparation

##### **Step 19.**

Turn on the vacuum pump and draw samples onto the membrane.

#### Membrane preparation

##### **Step 20.**

Wash each well with 400 µl of 5X SSC (2X).

#### Membrane preparation

##### **Step 21.**

Remove membrane and crosslink under UV light (200 mJ).

##### **Step 22.**

Store the membranes wrapped in plastic at 4°C.

#### **NOTES**

**Irina Agarkova** 05 Apr 2016

Note: Load the membranes so that when scanned, the second blot of the set may be placed under the first, i.e., load as the diagram in guidelines shows.