

Double SPRI for Second Generation Sequencing

Matthew Sullivan Lab

Abstract

Citation: Matthew Sullivan Lab Double SPRI for Second Generation Sequencing. **protocols.io**

dx.doi.org/10.17504/protocols.io.c53y8m

Published: 04 Jan 2016

Protocol

Step 1.

Add 65 μ l SPRI Ampure beads to 50 μ l of sample

Step 2.

Mix and incubate for 20 minutes

 DURATION

00:20:00

Step 3.

Separate on magnet (6 min.)

 DURATION

00:06:00

Step 4.

Transfer all supernatant (about 115 μ l) into a new well (off magnet)

Step 5.

Add 100 μ l of SPRI beads

Step 6.

Mix and incubate for 15 minutes

 DURATION

00:15:00

Step 7.

Separate on magnet (6min), discard supernatant

 DURATION

00:06:00

Step 8.

Wash beads with 60 μ l of 70% EtOH

Step 9.

Move well off magnet

Step 10.

Dry beads of EtOH (10min)

 DURATION

00:10:00

Step 11.

Add 40 µl EB

Step 12.

Mix and incubate for 3 minutes

 DURATION

00:03:00

Step 13.

Separate on magnet and transfer eluted product into destination plate

Warnings

Please note that depending on your application, it may or may not be necessary to eliminate high MW fragments assuming your shearing profile is relatively tight. And just eliminating the small material will help to increase yields. Make sure when you add the beads you mix well as this has dramatic effects on the yield.