



Sep 23, 2019

Taq PCR (Protocol for ZymoTaq™ PreMix)

[Alba Balletbó¹](#)¹Wageningen University[1](#) *Works for me* [dx.doi.org/10.17504/protocols.io.7cvhiw6](https://doi.org/10.17504/protocols.io.7cvhiw6) [Alba Balletbó](#)

GUIDELINES

ZymoTaq™ PreMix is supplied as a 2X concentrated “master mix” containing all the reagents needed to perform “hot start” PCR. The inclusion of a heat-activated, thermal-stable DNA polymerase reduces primer dimer and non-specific product formation that can occur when performing conventional PCR. This unique product is specifically designed for the amplification of bisulfite-treated DNA for methylation detection, realtime and quantitative PCR that are SYBR Green and probe based. The ZymoTaq™ PreMix yields specific amplicon formation with little or no byproducts. Simple and easy to use, just add water, primers, and template DNA to the ZymoTaq™ PreMix and then heat at 95°C for 10 minutes to initiate polymerization.

ZymoTaq™ DNA polymerase is a heat-activated, “hot start” polymerase that has 3'-terminal transferase activity. The addition of “A” overhangs to amplified DNA makes it ideal for use in TA-cloning.

Store ZymoTaq™ PreMix at -20°C for up to 12 months. Avoid repeated freeze/thawing of reagents. Prolonged storage is at -80°C

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
ZymoTaq™ PreMix - 50 rxns	E2003	Zymo Research
ZymoTaq™ PreMix - 200 rxns	E2004	Zymo Research

MATERIALS TEXT

Reagent	Volume	Final concentration
ZymoTaq™ PreMix	25 µL	1X
Forward Primer (10 µM)	Variable	0.3 to 1 µM
Reverse Primer (10 µM)	Variable	0.3 to 1 µM
Template	Variable	< 200 ng/50 µL
ddH2O	to 50 µL	
Total volume	50 µL	

- 1 Add all components in a 250 μ L tube making up to a 25 or 50 μ L reaction. If performing various PCR with different templates, a Master Mix is recommended to be done.

Reagent	Volume	Final Concentration
ZymoTaq™ PreMix	25 μ L	1X
Forward Primer (10 μ M)	Variable	0.3 to 1 μ M
Reverse Primer (10 μ M)	Variable	0.3 to 1 μ M
Template	Variable	< 200 ng/50 μ l
ddH ₂ O	to 50 μ L	
Total Volume	50 μ L	

- 2 Gently mix the PCR reactions and transfer the tubes to a thermocycler. Thermocycling conditions for a routine PCR:

Step	Temperature	Time	Cycles
Initial Denaturation		10 minutes	1
Denaturation	94-96°C	30 sec	30-40
Annealing	Variable	30-40 sec	
Extension	72°C	30 sec - 1 minute*	
Final Extension	72°C	7 minutes	1
Hold	4°C	> 4 minutes	1

* For \leq 1kb. Add an additional 15-30 seconds to the extension time for each kb > 1 kb. Make adjustments to the temperature and/or time if necessary.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited