

# Transformation of *Vibrio natriegens* Version 2

Carlos Helbig

## Abstract

This protocol describes how to transform chemically competent *Vibrio natriegens* cells.

The protocol was described and published by Weinstock et al., 2016

**Citation:** Carlos Helbig Transformation of *Vibrio natriegens*. **protocols.io**

dx.doi.org/10.17504/protocols.io.ptsdnne

**Published:** 30 Apr 2018

## Before start

You need:

- Chemically competent *V. natriegens* cells
- BHI + V<sub>2</sub>-Salts

### 10x V<sub>2</sub>-Salts:

| Chemical          | Concentration | Molecular weight | Amount for 1L |
|-------------------|---------------|------------------|---------------|
| NaCl              | 204mM         | 58.44g/mol       | 11.92176g     |
| KCl               | 4.2mM         | 74.55g/mol       | 0.31311g      |
| MgCl <sub>2</sub> | 23.14mM       | 203.3g/mol       | 4.704362g     |

### BHI:

For 1L medium 37g Brain heart infusion buillon

## Protocol

### Step 1.

Defrost your competent cells on ice

### Step 2.

Add your DNA template to the cells and mix gently

**Step 3.**

Incubate: 30min., on ice

**Step 4.**

Heatshock: 42°C, 45sec.

**Step 5.**

Incubate: 90sec., on ice

**Step 6.**

Add 1ml BHI + V<sub>2</sub>-Salts (pre-heated 30°C) to the cells

**Step 7.**

Incubate: 30°C, 1min.

**Step 8.**

Recover: 30°C, 120min., 200rpm

**Step 9.**

Centrifuge: 3000g, 4min.

**Step 10.**

Plate the pellet on warm agar plates (30°C)

**Step 11.**

Incubate: