



Oct 20, 2019

Pichia pastoris transformation through electroporation

[iGEM Dusseldorf¹](#)¹Heinrich-Heine Universität Düsseldorf[1](#) Works for me [dx.doi.org/10.17504/protocols.io.8heht3e](https://doi.org/10.17504/protocols.io.8heht3e)

iGEM Dusseldorf ⚡

MATERIALS TEXT

linearized vector DNA with electrocompetent Pichia cells

1 mL of 1:1 1 M sorbitol

YPD Liquid Medium or Plates (500 ml):

- 5 g yeast extract
- 10 g peptone
- 10 agar (if preparing plates)
- Fill to 450 mL with water

-> autoclave

YPDS + antibiotic Plates (1 liter):

- 10 g yeast extract
- 20 g peptone
- 182.2 g sorbitol
- Fill to 800 mL with water
- Mix until dissolved
- Transfer 400 mL two 1 L bottles
- pre-filled with 10 g BactoAgar
- Fill each bottle to 450 mL with water

-> autoclave

- Cool bottles to ~60°C, add 50 mL of 10x D
- Mix well, then aliquot in 4 bottles of 250 mL

Prior

- 1 Prior to transformation plasmid DNA containing the gene(s) of interest was linearized in restriction digests using an enzyme that only cuts once following the manufacturers manual !

Transformation

- 2 Mix approx. 150 ng of linearized vector DNA with electrocompetent Pichia cells
- 3 Transfer mixture to ice-cold 0.2 cm electroporation cuvette
- 4 Incubate cuvette for 2 min

- 5 Electroporate using a BioRad MicroPulser with a charging voltage of 1.5 kV
- 6 **Immediately** add 1 mL of 1:1 1 M sorbitol and YPD (v/v) to recover the cells
- 7 Transfer mixture to a plastic culture tube and incubated at 30°C for 1 hour at 225 rpm.
- 8 Centrifuge culture tube at 3.000 x g for 5 min to pellet the cells
- 9 Resuspend pellet in 200 µL YPD media
- 10 Plate out cells on YPDS plates containing desired antibiotic and incubated at 30°C for 3 days



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited