

# Sampling for Flow Cytometry (FCM) Version 2

Sarah Giuliani

## Abstract

*Collect and preserve samples for flow cytometry analysis. This protocol can be used for lab culture experiments or field samples.*

**Citation:** Sarah Giuliani Sampling for Flow Cytometry (FCM). **protocols.io**  
dx.doi.org/10.17504/protocols.io.f25bqg6

**Published:** 07 Oct 2016

## Materials

- ✓ P1000 micropipets and 1 ml filter tips by Contributed by users
- ✓ P200 micropipets and 200 µl filter tips by Contributed by users
- ✓ Labeled 1.2 ml cryo tubes by Contributed by users
- ✓ Repeater pipet and tips (that hold total volume of 100 µl and can dispense 2 µl) by Contributed by users
- ✓ Racks to hold cryo tubes by Contributed by users
- ✓ Cryo tube canes for dipping into liquid nitrogen by Contributed by users
- ✓ Labeled cryo boxes for -80°C storage by Contributed by users
- ✓ 25% Glutaraldehyde by Contributed by users
- ✓ Liquid Nitrogen: 1-3 L by Contributed by users
- ✓ artificial seawater salts by Contributed by users

## Protocol

### Step 1.

Advance preparation: for each sample to be collected, aliquot 900 µl of seawater (artificial seawater salts or complete sterile medium) into a labeled 1.2 ml cryovial.

### 📌 NOTES

**Maureen Coleman** 07 Oct 2016

*Adjust the dilution if necessary for your samples (this protocol uses a 10-fold dilution; if your sample is very concentrated and/or you are extremely sample-limited, you can dilute more; if your sample is dilute (<1E6 cells/ml), then you might want to dilute less or not at all).*

**Maureen Coleman** 07 Oct 2016

If your samples are not seawater-based, substitute an appropriate diluent of the same ionic strength (e.g. sterile medium, sterile freshwater salts).

## **Step 2.**

Collect sample: add 100 µl of sample into cryovial containing 900 µl of seawater. Keep cryovial cap clean and sterile.

## **Step 3.**

Preserve/fix sample: In the chemical fume hood, lay out cryovial caps (inside cap facing up) and use the repeater pipet to add 5 µl of glutaraldehyde (25%) to each cap. Cap each tube and mix well by inverting a few times. Incubate for 10 minutes in the dark.

## **NOTES**

**Maureen Coleman** 07 Oct 2016

Final concentration of glutaraldehyde is 0.125%.

## **Step 4.**

Flash-freeze tubes in liquid nitrogen and store at -80°C until analysis.