



Feb 06, 2020

## Synthesis of fluorinated neonicotinoids

PLOS One

Andrii Kyrylchuk<sup>1</sup>, Andriy Bezdudnyy<sup>1</sup>, Denys Klukovskyi<sup>1</sup>

<sup>1</sup>Institute of Organic Chemistry, National Academy of Sciences of Ukraine, Murmans'ka str. 5, Kyiv, 02660, Ukraine

1 Works for me dx.doi.org/10.17504/protocols.io.9h5h386

Andrii Kyrylchuk

### ABSTRACT

This is a general protocol for synthesis of fluorinated neonicotinoid analogues by interaction of amines with 2-chloro-5-(chloromethyl)pyridine.

### EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0227811>

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Mesquita RdS, Kyrylchuk A, Grafova I, Kliukovskyi D, Bezdudnyy A, Rozhenko A, Tadei WP, Leskelä M, Grafov A (2020) Synthesis, molecular docking studies, and larvicidal activity evaluation of new fluorinated neonicotinoids against *Anopheles darlingi* larvae. PLoS ONE 15(2): e0227811. doi: [10.1371/journal.pone.0227811](https://doi.org/10.1371/journal.pone.0227811)

### MATERIALS

NAME	CATALOG #	VENDOR
2-chloro-5-(chloromethyl)pyridine	516910	
anhydrous potassium carbonate	60109	
anhydrous acetonitrile	271004	
Dichloromethane	320269	Sigma Aldrich
Silica gel	645524	Sigma Aldrich
Ethyl acetate	319902	Sigma Aldrich
Hexane	208752	Sigma Aldrich

### MATERIALS TEXT

Substituted aniline, obtained according to literature sources.

- 1 Dissolve **2 mmol** of corresponding substituted aniline in **5 ml** of anhydrous acetonitrile in round-bottom flask, add **0.324 g (1.9 mmol)** of 2-chloro-5-(chloromethyl)pyridine and **0.828 g (6 mmol)** of anhydrous potassium carbonate.
- 2 Put a water-cooled backflow condenser on top of the flask and heat the flask under vigorous stirring until the solution starts to boil.

- 3 Continue refluxing for several hours. Monitor reaction progress via thin layer chromatography (TLC).
- 4 After the reaction is complete, filter off inorganic salts and rigorously wash the sediment with dichloromethane.
- 5 Combine all liquid fraction and remove the solvents by evaporation under reduced pressure.
- 6 Purify the residue by preparative (TLC) on SiO<sub>2</sub> (eluent EtOAc/hexane 1:2).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited