

# Top agarose

Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen

## Abstract

For use in "Isolation of cyanophages by plaque assays"

**Citation:** Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen Top agarose. **protocols.io**

dx.doi.org/10.17504/protocols.io.dqh5t5

**Published:** 02 Feb 2016

## Before start

Prepare 100 mL portions of 0.4 to 0.5% (w/v) of purified agar, agarose or low-melting point (LMP) agarose (i.e., Invitrogen #15517-022) in your media of choice. Although LMP agarose can be quite expensive, it is recommended for temperature sensitive samples and cells, since it solidifies at ca. 25°C. Purified agars, as well as low-melting point agars and agaroses are available for a range of lower temperatures (consult the following websites for more details: [www.sigmaaldrich.com](http://www.sigmaaldrich.com) and [www.invitrogen.com](http://www.invitrogen.com))

## Protocol

### Step 1.

Autoclave or microwave sterilize on the day of the assay.

### Step 2.

Dispense 2.5 to 3 mL into 13-x-100-mm glass disposable culture tubes (Fisher Scientific #1496127).

### Step 3.

Transfer tubes to a water bath or dry heat block set at the appropriate temperature.

#### ⊕ NOTES

**Amy Chan** 02 Sep 2015

(e.g., ca. 30 to 32°C for LMP agarose, ca. 40–42°C for purified agar or agarose)

### Step 4.

Cover tubes with foil or cap, allow for temperature to equilibrate.

### Step 5.

For each water sample, prepare triplicate tubes.

#### ⊕ NOTES

**Amy Chan** 02 Sep 2015

Control tubes containing cells are only used to monitor lawn growth.

**Amy Chan** 02 Sep 2015

For best results (smooth lump-free top agar), use freshly prepared top agar/agarose since repeated re-melting of solidified agar/agarose can give inferior results.