

RNA Samples Steps

Matthew Sullivan

Abstract

For [Transcriptomics During One-Step Growth Curves for Cellulophaga Phages protocol](#).

Citation: Matthew Sullivan RNA Samples Steps. **protocols.io**

dx.doi.org/10.17504/protocols.io.dev3e5

Published: 02 Feb 2016

Protocol

Step 1.

Use a serological pipet to sample 30 ml from each flask into a 50 ml falcon tube

Step 2.

Centrifuge for 5 min. at 4000 rpm, 4°C

 **DURATION**

00:05:00

Step 3.

Remove the supernatant

Step 4.

Resuspend the pellet in 1 ml [RNALater](#) and transfer to a 1.5 or 2 ml tube

Step 5.

Store samples at 4°C for up to 4 weeks