



Jul 02, 2019

IVF Bioscience Bovine Slaughterhouse Protocol

Jake Silcock¹¹IVF Bioscience

Working

dx.doi.org/10.17504/protocols.io.4rrgv56



Jake Silcock ⚡

Aspiration and IVM

1


Dish and Media Preparation

- 1.1 Preheat the Wash medium to 36°C
- 1.2 Prepare and equilibrate IVM wells at 38.8°C in 6% carbon dioxide in a humidified atmospheric air (21% oxygen) for 02:00:00 before use
- 1.3 Add $500\ \mu\text{l}$ of BO-IVM medium per well to 4-well plates without an Oil overlay and place in the incubator
- 1.4 Additionally, prepare one 35 millimetre dish with $3\ \text{ml}$ of BO-IVM per 4-well plate and place in the incubator without an Oil overlay







2

Aspiration

- 2.1 Set the heating plate to 34°C and connect the vacuum equipment
- 2.2 Warm the saline solution to 32°C
- 2.3 Take note of the arrival time and temperature of the ovaries
- 2.4 Wash the ovaries in the warmed saline and hold in a thermo-container of warmed saline
- 2.5 Add $140\ \mu\text{l}$ of heparin to the 50 ml tube prior to aspiration

- 2.6 Aspirate all follicles between 2-15 millimetres using an 18-guage needle connected to the pump system. The opening of the needle should face downwards
- 2.7 Aspirate no more than  20 ml of follicular fluid per tube and note the number of ovaries aspirated on the lid of the tube

3 Collection and Maturation

- 3.1 Use 9 cm Petri dishes marked with a square grid pattern
- 3.2 Prepare one dish per tube of follicular fluid. Mark the dish with the corresponding tube number.
- 3.3 Add  7 ml of Wash to each dish
- 3.4 Prepare three 35 millimetre dishes with  2.5 ml of Wash in each. Mark these dishes as 1, 2 and 3
- 3.5 Transfer the oocyte pellet from the bottom of the 50 ml tube into the larger petri dish using a plastic transfer pipette
- 3.6 Search through the dish systematically and transfer all oocytes to the smaller dish marked 1
- 3.7 Repeat the transferring of the oocyte pellet and search, twice for each tube
- 3.8 Once all oocytes have been collected, wash through dish 2 and then 3. Count the number of collected oocytes
- 3.9 Rinse once through the 35 millimetre dish containing BO-IVM
- 3.10 Once rinsed, transfer the oocytes to the 4-well plate containing BO-IVM. Do not transfer more than 45 oocytes per  500 µl well.
- 3.11 Incubate at  38.8 °C and 6% carbon dioxide in humidified atmospheric air (21% oxygen) for between  21:00:00 and  24:00:00 .



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited