

Culturing Symbiodinium

Tingting Xiang, Arthur Grossman

Abstract

we generated five clonal, axenic strains of Symbiodinium. These strains were assigned to clades A (two strains), B, E, and F based on their chloroplast 23S rDNA

sequences. Growth studies in liquid cultures showed that the clade B strain and one of the clade A strains were able to grow photoautotrophically (in light with

no fixed carbon), mixotrophically (in light with fixed carbon), or heterotrophically (in dark with fixed carbon). Here this protocol describes the steps of culturing SSB01 in different types of medium as an example.

Citation: Tingting Xiang, Arthur Grossman Culturing Symbiodinium. protocols.io

dx.doi.org/10.17504/protocols.io.qg8dtzw

Published: 23 Jul 2018

Materials

- \checkmark CORAL PRO SALT R_99344 by Contributed by users
- DAIGOS IMK MEDIUM FOR MARINE MICROALGAE 398-01333 by Wako
- Glucose G8270 by Sigma Aldrich
- Marine Broth 2216 279110 by BD Biosciences
- Noble agar A5431 by Sigma Aldrich
- Enzymatic Digest Casein Hydrolysate AAJ12855Q1 by Contributed by users

Protocol

Prepare liquid media

Step 1.

Artificial seawater (ASW) is the base for most of the media, except Marine Broth.

Making ASW

Step 2.

Dissolve 37.1 g Coral Pro salt into 1L pure water; check the salinity of the dissolved seawater and make sure it is between 32-34.

₽ NOTES

This step of dissolving sea salts may take a while, stir overnight is recommended.

Making 10x Daigo stock solution for the media

Step 3.

Dissolve 2.56g Daigo's IMK Medium powder into 1L ASW;

₽ NOTES

Daigo's IMK Medium is the base of IMK medium.

Making 10x Daigo's IMK Medium stock solution

Step 4.

Filter sterile the 10X Daigo's IMK Medium solution through a 0.22 µm filter;

P NOTES

Here is the recipe of Daigo' IMK Medium, which contains nutrients such as vitamins. Therefore, filter sterilization is highly recommended instead of autoclave sterilization.

Published: 23 Jul 2018

[IV]細胞培养用培地

微細蒸類培養用

海座微細藻類培養用

ダイゴIMK培地

Daigo's IMK Medium for Marine Microalgae

【和光コード・包装】398-01333 (100L用×10) 392-01331 (1000 L用)

概要

本培地は、多大な機能藻類をできるだけ多く確実に効率よく 培養可能とする培地を簡便に調製できるよう、「株式会社 海洋パイオテクノロジー研究所」が開発した海産微細薬類用

の培地です。 海水に溶かすだけで、培地の調製が可能なように、全ての要 素物質が混合されています。

環境中の微細藻類の分離から大量培養(性)まで広覧に利用可 能で、海産微細藻類培養、生理学的研究、水産業における種 苗生産用の飼料藻類の培養等に便利な培地として設計されて

います。 (池)建漢類の高密度培養には、別途建酸ナトリウム(0.2~LmM)を添加

組成・成分 (mg/L)

| NaNO ₃ | 200 | CoSO ₄ · 7H ₂ O | 0.014 |
|---------------------------------------|-------|--|--------|
| Na ₂ HPO ₄ | 1.4 | Na ₂ MoO ₄ • 2H ₂ O | 0.0073 |
| K ₂ HPO ₄ | 5 | CuSO ₄ · 5H ₂ O | 0.0025 |
| NH ₄ C1 | 2.68 | H ₂ SeO ₃ | 0.0017 |
| Fe-EDTA | 5, 2 | Thiamin-HCl | 0. 2 |
| Mn-EDTA | 0.332 | Biotin | 0,0015 |
| Naz-EDTA | 37. 2 | Vitamin B ₁₂ | 0,0015 |
| ZnSO ₄ • 7H ₂ O | 0.023 | MmCl ₂ · 4H ₂ O | 0.18 |

使用法

海水100L (1000L) に本品25.2g (252g) を撹拌しながら加え て溶解する。

- 奈: pは耐きを必要とする場合は、塩酸又は水酸化ナトリウムを用いる。pはボアルカリ側に僅くと自色沈酸を生することがある。
 奈: ダイゴIMK培地を用いて大量培養や高密度培養を行う場合は、11に対する粉末原加量を増やすことが可能である。
- ※: IMK培地を人工海水SPに溶解するとき、必ず人工海水SPを溶解 した後にIMK培地を溶解してください。
- 1) 滅菌について
- 減蓄をしていない培地は、培地調製中に海水または空中の機組藻類が混入し、治培所においても微組藻類が発生することがある。
- ダイゴIM K培地には、ビタミン類が調合されているので、培地を被菌する際には、培養対象のビタミン要求性を考慮して、メンブランフィルターを用いてろ追ば菌する。減菌後は4℃で保管し、できるだけ早く使用すること。
- 2)カンテン培地の調製法

カンナン暗地の調発法 ダイゴ人工権水野(管本) 500mLにダイゴIMK培地202mgを加え海 解し、減菌後地地を50で程度に保つ。別に、素留水200mLにカンテン 15gを加え121で19分間高圧蒸気破離し、提井の後、幸風で放置し で初で程度になったら、先に破菌した増地とよく混合して、シャー レに分性する。固化後、希腊所に保管する。

注意事項

※本品は『医薬用外毒物』で

本品は吸湿性が強いので開封後は密封して保管すること。 室温保存。

32

Making liquid IMK medium

Step 5.

Transfer 90 mL ASW into a 250 mL PYREX Glass Erlenmeyer Flask, and cover with aluminum foil;

Making liquid IMK medium

Step 6.

Autoclave at 121 °C for 30 minutes;

Making liquid IMK medium

Step 7.

Let it cool to room temperature.

Making liquid IMK medium

Step 8.

Bring the flask from step 7 and 10X Daigo's IMK Medium stock solution from step 4 to a sterile hood;

Making liquid IMK medium

Step 9.

Add 10 mL Daigo's IMK Medium stock solution to the flask; and that makes 100 mL IMK.

Making liquid IMK medium supplimented with casein hydroplysate (IMK+cas)

Step 10.

Dissolve 0.4 g casein hydroplysate into 90 mL ASW, and transfer to a clean 250 mL PYREX Glass Erlenmeyer Flask, and cover with aluminum foil;

₽ NOTES

Please note to use enzymatic digest casein hydrolysate, the hydrochloric acid hydrolysate of casein may not dissolve well in ASW.

Making liquid IMK medium supplimented with casein hydroplysate (IMK+cas)

Step 11.

Autoclave at 121 °C for 30 minutes;

Making liquid IMK medium supplimented with casein hydroplysate (IMK+cas)

Step 12.

Let it cool to room temperature.

Making liquid IMK medium supplimented with casein hydroplysate (IMK+cas)

Step 13.

Bring the flask from step 12 and 10X Daigo's IMK Medium stock solution from step 4 to a sterile hood;

Making liquid IMK medium supplimented with casein hydroplysate (IMK+cas

Step 14.

Add 10 mL Daigo's IMK Medium stock solution to the flask; and that makes 100 mL IMK+cas medium.

Making liquid IMK medium supplemented with glucose (IMK+Glc)

Step 15.

Dissolve 0.5 g glucose into 90 mL ASW, and cover with aluminum foil;

Making liquid IMK medium supplemented with glucose (IMK+Glc)

Step 16.

Autoclave at 121 °C for 30 minutes;

Making liquid IMK medium supplemented with glucose (IMK+Glc)

Step 17.

Let it cool to room temperature.

Making liquid IMK medium supplemented with glucose (IMK+Glc)

Step 18.

Bring the flask from step 17 and 10X Daigo's IMK Medium stock solution from step 4 to a sterile hood;

Making liquid IMK medium supplemented with glucose (IMK+Glc)

Step 19.

Add 10 mL Daigo stock solution to the flask; and that makes 100 mL IMK+Glc medium.

Making liquid Marine Broth (MB) medium

Step 20.

Dissolve 37.4 g MB powder into 1L pure water;

Making liquid Marine Broth (MB) medium

Step 21.

Transfer the liquid to a new 1L glass container;

Making liquid Marine Broth (MB) medium

Step 22.

Autoclave at 121 °C for 30 minutes; and let it cool down to room temperature.

P NOTES

There might be some precipitation after autoclave, filter through a 0. $45\mu m$ filter is recommended to get rid of the precipitation.


Inoculate the Symbiodinium culture

Step 23.

Pick single colonies from agar plates which contain the Symbiodinium colonies;

P NOTES

Image showing Symbiodinium colonies on a solid plate.



Inoculate the Symbiodinium culture

Step 24.

Grow the cultures at 27°C under 10 μmol photons $m^{\text{--}2}\,s^{\text{--}1}$ white light.

₱ NOTES

